Observation

Stéréomicroscopes

Observation

Stereomicroscopes


Réf :
571 340

Français – p 1

English – p 5

Version : 3107

Loupe binoculaire autonome
Binocular stereomicroscope

Légende :

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
0825 563 563 *

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Caption :

- | | |
|---|---------------------------------------|
| A Wide field eyepieces with eyeshades | E Base with built-in batteries |
| B Dioptic adjustment | F Binocular head |
| C Objectives knob | G Macrometric adjustment |
| D LED incident and/or transmitted illumination | H Illumination switches |
| | I Brightness adjustment knob |

1. Presentation

This binocular stereomicroscope is meant to observe any biological or geological material, which cannot be observed with a microscope due to its thickness and/or its opacity.

This is a stand-alone device with its rechargeable batteries (supplied with). It is a model with both incident and transmitted illumination.

2. Specifications

| | |
|--------------------------------|---|
| Head | Binocular, 45° inclined, 360° rotatable |
| Eyepieces | Wild field x10 with eyeshades |
| Interpupillary distance | 55 - 75 mm |
| Dioptric adjustment | On the 2 tubes |
| Objectives | x2 and x4 |
| Magnification | x20 and x40 |
| Field of view | 11 mm / x20 ; 5,5 mm / x40 |
| Working distance | 100 mm |
| Stage | 1 white/black and 1 transparent |
| Illumination | Both incident and transmitted LED |
| Power supply | Rechargeable batteries (supplied with) and 4,5 V / 600 mA transformer |
| Dust cover | Yes |

3. Unpacking and assembling

3.1. Unpacking the stereomicroscope

Carefully open the styrofoam packaging and check that all components have been taken out the packaging. Do not throw out the packaging before identifying all the components.

3.2. Assembling the stereomicroscope and first use

When manipulating the components, and particularly optical parts, take care to not touch the lenses with your bare fingers or hands. Any fingerprints or grease spots will reduce the image quality.

After unpacking the support, install it on a stable working surface and remove the various protections items (do not touch the surface of the lens). Fix the head and install the eyeshades on the eyepieces.

It exists 2 ways to supply the stereomicroscope:

- stand-alone with the supplied rechargeable batteries,
- with the transformer.

First use: Before connecting the instrument and switching it on, make sure that the operating voltage of the stereomicroscope is the same as the power supply voltage.

Batteries specifications:

Loading time = 4 hours; Stand-alone time = 4 hours.

3.3. Alignment and working

3.3.1. Interpupillary distance

Move the two eyepiece tubes until only one circular field is visible through the two eyepieces. If two circles appear, the interpupillary distance is too large, and if 2 overlapping circles appear, the interpupillary distance is too small.

3.3.2. Focusing

Turn the objectives selection's knob to the biggest magnification. Bring the specimen into focus until the image is clear and sharp with the macrometric knob.

Turn the objectives selection's knob to the smallest magnification. Adjust the dioptic knob of the right eyepiece until the image of the right eyepiece is clear and sharp.

Repeat the procedure for the left eyepiece. Then check the image focus over the entire magnification range. It should then be **perfectly parfocal**.

3.3.3. Illumination

For your comfort, this stereomicroscope is supplied with separate incident and transmitted illumination. Adjust the brightness with the illumination knob.

4. Cleaning and maintenance

Protection against dust: Cover the instrument with the protective cover provided if it remains unused for a long period. Never leave the eyepiece tube open and unprotected. Either leave the eyepiece in the tube (which is recommended if the microscope is used very frequently) or cover it with wrapping paper or a protective cap. It is recommended that eyepieces and other optical accessories should be kept in a dry box to protect them from dust and humidity when they are not in use.

Protection against humidity: The instrument must be kept at a distance from every water source, pipe and sink. Humidity in the room in which the instrument is to be used must be as low as possible (the relative humidity must remain less than 70%). It is recommended that all optical accessories should be kept in a dry box when they are not being used. It is strongly recommended that a dehumidifier and air conditioning should be used continuously if the environment is too humid.

Cleaning: If there is any dust present on the optical surface, try to remove it with a fan or compressed air.

There are two recommended possible methods for removing fingerprints, grease or dust stains that cannot be removed by blowing with air :

- Blow gently onto the glass surface and then wipe with a clean cloth, paper for lenses or a cotton swab. Note that small cotton fibres may remain on the surface of the lens if a cotton swab is used.
- Use a cotton swab to clean the surface of the lens thoroughly. Do not use aggressive solvents.

The shelf life of the LED is of 5,000 hours. In case of non-function of the lighting and after checking that the batteries are charged, send the stereomicroscope back to Jeulin.

5. Moving the stereomicroscope

The microscope must not be moved if it can be avoided. If it must be moved, the user must check that the eyepieces are well fixed into the eyepiece tubes before moving the instrument.

6. After-sales service

The guarantee is 2 years. The equipment must be returned to our workshops.

If you need any repairs, adjustments or spare parts, please contact:

JEULIN - Technical support
468 rue Jacques-Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries