Observation

Microscopie - Accessoires

Microscopy

Microscopy - Accessories

Réf :
574 078

Français – p 1

English – p 4

Version : 5109

Platine polarisante
Polarizing stage

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27019 EVREUX CEDEX France

0 825 563 563*

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The polarizing stage is an easy-to-use product allowing to **adapt your biological microscope or your stereomicroscope (with incident and transmitted illumination)** to geology applications.

This is a unique device offering an excellent optical and mechanical quality.

**Don't comply with Naja microscopes and Irio with double layers mechanical stage.*

1 Aim and description

Together with a biological microscope, the polarizing stage allow to observe cristalline structures in polarizing light like starch, aspirin but mainly in geology, rock slides.

1.1 Specifications

This polarizing stage consists of :

- 1 retractable analyzer,
- 1 graduated disc over 180° by 5° allowing to place the rock slide and change the extinction angle of minerals,
- 1 polarizer.

Dimension of the stage: 98 x 40 mm.

Diameter of the graduated disc: 70 mm.

1.2 Application

Use an 80x magnification max. (as example 15x eyepiece together with 4x objective) only to get the larger field of view.

Place this device on the stage of your microscope or stereomicroscope (with incident and transmitted illumination) and maintain it with the clips or the slide-guide (depends on the product you use).

Select the lowest magnification. Focus the device turning the graduated disc and try to keep the minerals in the field of view.

Remove the analyzer (upper part of the device) and put the rock slide on the center of the graduated disc.

Turn the graduated disc: the central mineral should stay in the middle of the field of view. If necessary, move the device on the stage, turn again the graduated disc to get the minerals centered in the field of view. Try again with an higher magnification: the polarizing stage is now centered.


Remove the rock slide.

Put back the analyzer and turn it to get no light towards the system. The analyzer and the polarizer are crossed.

Your microscope or stereomicroscope is ready to use.

1.3 Aim

Light Passing Through Crossed Polarizers

The concept illustrated in the above drawing explores the interaction of light waves with polarizers. In this example, the incident light is vibrating perpendicular to the direction of propagation in an equal distribution of all planes before encountering the first polarizer. The polarizers illustrated above are oriented in a single direction. Only the incident light that is vibrating in the same plane is absorbed, while light vibrating at right angles to the plane is passed through the first polarizing filter. Polarizer 1 (analyzer) is oriented vertically to the incident beam so it will pass only the waves that are vertical in the incident beam. The wave passing through polarizer 1 (analyzer) is subsequently blocked by polarizer 2 because the second polarizer is oriented horizontally with respect to the light wave. The concept of using two polarizers oriented at right angles with respect to each other is commonly termed ***crossed polarization*** and is fundamental to the practice of polarized light microscopy.

When placing the rock slide between the 2 polarizers, the polarizing plane turn because of the slide which refract the polarizing light. But the polarizing plane is not perpendicular to the polarizer 2 anymore. The slide is creating interferences in a non-homogeneous environment because the beam go through minerals with different speed.

1.4 Pratique

When polarizing stage ready to use, place the rock slide in the center of the graduated disc. Take care to not move the analyzer and observe the colors of minerals.

When turning the graduated disc, you can observe the colors of minerals change four times. The angle when the color change depends on the mineral.

2 After sales department

The device is under a 2 years guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux