Chromatographie

Chromatography

**Ref :
591 029**

Français – p 1

English – p 5

Version : 6011

Kit Chromatographie

Chromatography Kit

1 Objectif

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1 Purpose

This chromatography kit enables you to perform chromatographies on paper, in particular the separation of photosynthetic pigments. Chromatography on analytic filter paper strips is less costly and easier to reproduce than on TLC.

Other applications: chromatography of sugars, amino acids, polyfunctional compounds.

Summary of the theory

Liquid-solid chromatography is based on the separation of solutes between the stationary adsorbant and the mobile liquid. Each solute is subject to a retention force (by adsorption) and a transport force by the mobile phase. The equilibrium between these two forces results in a differential migration of the solutes in the sample being analyzed in relation to their molecular characteristics, the size and polarity of the molecule, which causes them to separate. Adsorption is the fixation of the dissolved molecules by the chromatography support (paper or silica gel). This fixation results from the establishment of weak chemical bonds; e.g., dipolar ionic, dipolar-dipolar, or Van der Waals bonds, between the adsorbant and the molecule adsorbed.

2 Contents

6 vessels, 6 stoppers with hooks, 1 set of 36 strips of analytic filter paper (24 x 285 mm).

- Vessel: Stable base and height suitable for the separation of chlorophylls, xanthophylls and carotenoids. Without a spout, and ungraduated to optimize saturation and readability.
- Stopper: The design of the stopper makes it possible to fasten the chromatography paper and saturate the vessel with vapors. Its antishock rim makes it easier to grip and protects the vessel from breaking in the event of a spill, and it also prevents the dispersion of the solvent, which is potentially toxic.

3 Examples of applications


The study of chlorophyll pigments: separation of the pigments by chromatography

The first chromatography was conducted in 1906 by the Russian botanist M. TSWETT. He used it to separate the pigments in a spinach leaf. M. Tswett observed the separation of the plant pigments composing whole chlorophyll when he filtered a chlorophyll solution in petroleum ether on a calcium carbonate column. This resulted in the formation of green and yellow colored zones. M. Tswett concluded: "Exactly like the light rays in the spectrum, the different components of a mixture of pigments separate out on the calcium carbonate column according to a law and they can be analyzed qualitatively and quantitatively." This invention of chromatography remained unknown until 1931, when it was successfully used again. Different supports are used: paper, metal plates coated with silica gel, etc. The location of the pigments is different for each type of support, paper or silica gel.

3.1 Instructions

3.1.1 Extraction of foliar pigments:

- Finely cut 5 healthy, clean ivy leaves in a mortar containing a pinch of sand and 20 ml of denatured ethanol.
- Mash with a pestle until a thick, green liquid is obtained,
- Filter the liquid thus obtained to recover an alcoholic solution of foliar pigments.

3.1.2 Preparation of the solvent:

To make 100 ml of solvent, mix 85 ml of petroleum ether, 10 ml of acetone, and 5 ml of cyclohexane.

3.1.3 Separation of the pigments by chromatography:

Preparation of the chromatography vessel:

Add extraction solvent to the chromatography vessel to a depth of 1 cm, then stopper the vessel to saturate it with the solvent vapor.

Preparation of the support:

- Draw a thin line in pencil 2 cm from one end of the strip.
- In the center of the line, deposit 6 to 10 drops of the alcoholic pigment solution, be sure to deposit small drops and let each drop dry before depositing the next one. The diameter of the final deposit must not exceed 8 mm ①.
- Carefully hang the strip on the hook on the vessel stopper.
- Introduce the strip into the vessel (the solvent must not cover the deposit).
- Let the solvent migrate and wait until the level of the solvent is within 10 mm of the upper end of the chromatography paper strip ③.
- Stop the process, remove the paper strip from the vessel and let it dry.

Note: The supports supplied with the chromatography kit vessels provide space to make two deposits on the same strip so that they can be compared.

①

②

③

3.2 Expected results:

This process will reveal four types of pigments in the following order: chlorophylls b and a, xanthophylls and carotenes ③.

3.3 Variations

This technique may be used for other biological materials such as: spinach, chard, algae, euglenophyta and cyanobacteria (blue-green algae), any plant with green leaves.

3.4 Tricks and tips

- Let the deposits dry well.
- Once the deposit is made, put the strip in the vessel quickly.
- Ideally, let the migration take place in a **dark place**, as the pigments are photosensitive.
- Take the proper precautions when handling solvents (masks, gloves, fume hoods, etc.).

Accessories	
Chemical products for the solvent	
Petroleum ether 900 ml	Ident. no. 102 008
Acetone (propanone) 1000 ml	Ident. no. 102 021
Cyclohexane 1L	Ident. no. 102 060
Analytic filter paper	
Strips 24 x 285 mm. Set of 36	Ident. no. 703 077
Whole sheets 250 x 285 mm. Set of 25	Ident. no. 703 076
Elkay® liquipettes	
Polyethylene spatula. L: 130 mm, Vol : 0.5 ml, 40 drops per ml, Vol. 1 drop: 25 µl.	Ident. no. 723 042
Fine tips. Set of 25.	

4 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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Service gratuit * :
+ 33 (0)2 32 29 40 50

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Rubrique FAQ

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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