Appareil de laboratoire

Congélation

Laboratory apparatus

Freezing

Réf :
701 075

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Mini-freezer

Français – p 1

English – p 19

Version : 6009

1 Généralités

1.1 But de l'appareil.

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1 General information

1.1 Purpose of the apparatus.

The mini-freezer is an apparatus designed to cool a liquid contained in a small 4 ml chamber. It was developed especially for studying changes of state of pure substances such as water or cyclohexane, or mixtures (water-salt, water-ethanol).

Equipped with an ergonomic thermometer support, this apparatus enables one to follow the progression of the temperature in the medium being studied during the change of state with the aid of a liquid or digital thermometer, or with the aid of a computer-assisted experimentation data acquisition chain.

1.2 Description

1.2.1 Views of the apparatus



Figure 1

- 1 – Inclinal thermometer support
- 2 – Function indicator light.
- 3 – Protective cover for the refrigeration device
- 4 – Body of the apparatus.
- 5 – Heat reservoir tub.
- 6 – Sliding, detachable thermometer support adapter.
- 7 – Refrigeration chamber.

- 8 – Cap for the heat reservoir fill hole.
- 9 – Polarity guide for the power supply to the apparatus.
- 10 – Handle for picking up and moving the apparatus.
- 11 – Depression permitting various units to be stacked.
- 12 – 4 mm Ø safety sockets for the power supply to the apparatus.

1.2.2 Refrigeration device

This apparatus is equipped with a refrigeration device that, when supplied with electricity, absorbs the heat from the refrigeration chamber (*Figure 1-7*) and disperses it in a large capacity heat reservoir (*Figure 1-5*) containing water that gradually heats up (*Figure 2*).

1.3 Technical specifications

Power drawn	40 W
Minimum temperature attained	-17 °C* *(For a heat reservoir in which the temperature ≤ 14 °C at the start of the experiment)
Power supply	12 V DC / 4 A
Connection	4 mm $\varnothing$ safety sockets. This apparatus is protected against polarity inversions.
Overheating protection	Resettable circuit breaker after freezing and power cutoff.
Function indicator	Green light indicating that the apparatus is working.
Dimensions (L×w×h) of the apparatus Weight (empty)	167×180×160 mm 600 g (ca.)
Refrigeration chamber: Construction material Capacity	Cast aluminium 4 mL
Heat reservoir: Construction material Capacity	Antishock polypropylene. 1.5 L
Agitator of refrigeration chamber	The mini-freezer is supplied with a manual agitator made of antishock polypropylene.
Thermometer support: Construction material Adapters	Antishock polypropylene. For all thermometers with $\varnothing 3$, $\varnothing 3.5$, $\varnothing 6$ mm

2 Startup

2.1 Use precautions

2.1.1 Power supply specifications

The mini-freezer must not be operated at a voltage greater than 12 V.

A voltage greater than 12 V may destroy the refrigeration device.

This apparatus is also protected against electrical connections that do not correspond to the polarity indicated on the apparatus (*Figure 1-9*).

2.1.2 Heat reservoir

The heat absorbed from the refrigeration chamber must be removed. **It is therefore essential to make sure that the heat reservoir is filled with sufficient water before using the apparatus.** If it is not, the apparatus may overheat and trip the emergency circuit breaker.

Do not put any liquid other than water in the heat reservoir.

2.1.3 Refrigeration chamber

The refrigeration chamber is made out of aluminum. For this reason, **it must not come into contact with acids, bases, or corrosives, or it may be permanently damaged.**

2.1.4 Spills

If water is spilled on the apparatus, wipe it off with an absorbant tissue or a sponge.

Never immerse the mini-freezer in water.

2.1.5 Optimum operating conditions

The operating principle of the apparatus is based on the transfer of the heat absorbed in the refrigeration chamber to the heat reservoir. The output of the apparatus (and thus its efficiency) diminishes as the water in the heat reservoir heats up. It is therefore suggested to fill the heat reservoir of the apparatus with cold tap water shortly before initiating a series of experiments. Using water that is too warm may prolong the time required to complete the experiment and increase the minimum temperature that can be obtained in the chamber.

2.1.6 Moving the apparatus

To move or pick up the apparatus, hold it by the handle (*Figure 1-10*).

Never pick up the apparatus by its body (*Figure 1-4*) or by the thermometer holder (*Figure 1-1*).

2.2 Setting up the apparatus

2.2.1 Filling the heat reservoir

Open the fill hole (*Figure 1-8*) of the heat reservoir of the mini-freezer as indicated below.


Figure 3


Figure 4

With the tip of your finger, press lightly on the side of the sprocket (*Figure 3*) located on the cap in order to loosen the latter. Slide the cap by pushing the sprocket until it goes in the slot (*Figure 4*).

Fill the reservoir by holding the fill hole under a thin stream of water from a faucet. Fill until the surface of the water reaches the juncture between the body of the apparatus and the heat reservoir tub (the division between the beige and brown colors).

Close the fill hole by pushing on the side of the sprocket located on the cap in order to slide the latter back into its original position.

2.2.2 Inserting the thermometer

The mini-freezer is equipped with a thermometer holder (*Figure 1-1*) to hold a liquid or digital thermometer, or a temperature probe in the freezing chamber. The thermometer holder also makes it possible to center the thermometer exactly in the freezing chamber without having it touch the walls and giving a false temperature reading. In the storage position, the thermometer holder is pushed down into a depression located on the protective cover of the refrigeration device (*Figure 1-3*).

To insert the thermometer, turn the thermometer holder to its upright position by pivoting it about its axis (*Figure 5*).


Figure 5

Turn the thermometer adapter a quarter of a turn (*Figure 6*).


Schéma 6

Slide the thermometer adapter along the post of the holder until the desired height is obtained (*Figure 7*).


Figure 7

The apparatus comes with three interchangeable thermometer adapters with the following diameters: $\varnothing$ 3, $\varnothing$ 3.5 or $\varnothing$ 6 mm. To change an adapter, remove the old adapter from the holder by sliding it off the post of the holder and then insert the adapter-support with the desired diameter (*Figures 8, 9, and 10*).


Figure 8
 3 mm $\varnothing$ adapter


Figure 9
 3.5 mm $\varnothing$ adapter


Figure 10
 6 mm $\varnothing$ adapter

There is a corresponding adapter for each type of thermometer or temperature probe. The table below lists the main types of thermometers available and the corresponding adapter (the list is not exhaustive).

Adapter	Thermometer or probe
3 mm Ø adapter	- VTT rapid response thermometric sensor
3.5 mm Ø adapter	- "Checktemp" digital thermometers
6 mm Ø adapter	- Liquid thermometers - Initio® -20/+120 °C thermometer - Initio® pH meter thermometer - VTT® -20/+120 °C thermometric sensor - Pt1000 ESAO® thermometric sensor

Clip the thermometer onto the adapter as shown below (*Figure 11*).


Figure 11

Figure 12

Do not attach the thermometer to the adapter as shown in Figure 12. Doing so may result in injury to the user or damage to the thermometer.

Make sure that the thermometer is immersed in the refrigeration chamber but not touching the bottom (*Photo 13*).

Photo 13

2.2.3 Power supply

Connect the mini-freezer to a 12 volt DC transformer via cables equipped with 4 mm diameter banana safety plugs (*Photo 14*). The power supply must furnish at least 4 amps of current. We highly recommend using the 6/12 volt Evolution® F6F12/5A DC transformer (Item no : 281 083) or a power source with the same specifications.

Make sure to observe the polarity indicated on the protective cover of the apparatus (*Figure 1-9*).

Turn on the power. The green indicator light (*Figure 1-2*) should come on, indicating that the apparatus is functioning.

If the indicator light does not come on, refer to the instructions in chapter 6.

Photo 14

2.3 Agitating the medium

In order to obtain good results while studying the changes of state by measuring the temperature as a function of time, it is essential to homogenize the temperature of the medium by agitation.

For this purpose, the mini-freezer is equipped with an agitator for homogenizing the temperature of the medium being studied. The agitator will not come into contact with the thermometer and using it avoids splashing and spilling the liquid contained in the refrigeration chamber.


Figure 15

Place the agitator in the refrigeration chamber so that the fins of the agitator are immersed in the liquid and so that the crown of the agitator rests on the rim of the chamber (*Figure 15*).

Grasp the crown of the agitator between your thumb and index finger and work the agitator by rotating it back and forth (*Figure 15*).

2.3.1 Filling the refrigeration chamber

To fill the refrigeration chamber:

- With a pipette or dropper, add the liquid to be studied into the refrigeration chamber.
- Immerse the tip of the thermometer in the liquid **without letting it touch the bottom of the refrigeration chamber** and place the agitator in the chamber.
- Adjust the level of the liquid until it comes to within 1 to 2 millimeters of the rim of the refrigeration chamber.
- If necessary, remove any excess fluid with the aid of a pipette or dropper.

To change the liquid:

- Remove the thermometer and the agitator from the refrigeration chamber, rinse them and wipe them dry.
- Empty the refrigeration chamber by aspirating its contents with a pipette or dropper.
- If necessary, rinse with distilled water, using a pipette or dropper.
- Then wipe up the residues by introducing a piece of absorbant tissue into the refrigeration chamber.
- Now you can fill the chamber with another liquid.

Using a dropper or a pipette avoids accidental spills of large quantities of liquid on the cover of the apparatus.

2.3.2 Prolonged functioning

After the solidification of the liquid in the refrigeration chamber occurs, the temperature will fall until it attains a minimum. At this point, turn the apparatus off. If it is allowed to run for a long time, the heat reservoir will heat up and the apparatus will no longer be able to evacuate the heat absorbed from the refrigeration chamber. The overheating protection mechanism may then be triggered, shutting off the power to the apparatus. The indicator light will then go out (*Figure 1-2*).

To reset the overheating protection mechanism, refer to paragraph 6.2.2.

3 Manipulations and experiments

3.1 General information

The mini-freezer is an apparatus designed for studying the changes of state of pure liquids and mixtures. With it, one can chill a small volume of the liquid being studied and study the progression of the temperature of the liquid as a function of time.

It is thus possible to study the behavior of the liquid as it chills and observe the shape of the curve $T = f(t)$.

The presence or absence of a temperature plateau during the solidification process allows one to address the idea of a pure substance or a mixture and to determine the freezing temperature of the liquid in either case.

Various types of liquids can be studied:


Pure substances	Demineralized water Cyclohexane
Mixtures	Demineralized water + salt (e.g., NaCl) Water + miscible liquid (e.g., ethanol). Automobile coolant.

3.2 Measurement precision

The precision of the results obtained depends mainly upon the type of thermometer used. In fact, the small volume of the study sample may be perturbed substantially by the measuring apparatus used. For this reason, the value actually measured may differ by several degrees from the theoretical value if the wrong thermometer is used.

It is therefore recommended to use thermometers or temperature probes that only perturb the system being studied very slightly (see the following table).

Manipulation	Suggested measuring apparatus
Manual measurements	"Checktemp" digital thermometers.
Computer-assisted measurements	VTT [®] console equipped with a rapid response thermometric sensor. ESAO [®] Primo console equipped with a rapid response thermometric sensor.

(Photo 16)

3.3 Tips and tricks

3.3.1 Observation of the solid formed

It may be interesting to observe the solid formed by freezing. To do so, turn off the power supply to the mini-freezer. Wait a few minutes until the refrigeration chamber warms up in contact with the ambient air (as indicated by the melting of the frost formed on the outer wall of the chamber during the experiment). Turn the temperature probe fixed in the solid mass, which should make it come loose from the wall of the refrigeration chamber. Carefully pull the temperature probe up to remove the solid mass from the chamber and observe it (Photo 16).

3.3.2 Reducing the power of the apparatus

When studying the change of state of certain substances having a low thermal capacity and a low enthalpy of fusion, it may be necessary to reduce the power of the mini-freezer.

With such substances, the duration of the experiment may otherwise be too short to allow a good observation of the phenomena taking place. It may also be difficult to take measurements. To remedy this problem and increase the duration of the experiment, supply the mini-freezer with 6V/5A of power.

As a power source, we recommend the Evolution[®] F6F12 5A transformer, which can supply 6 or 12 volts, depending on the refrigeration power required.

3.3.3 Optimizing the power of the apparatus

In order to obtain temperatures lower than -13 °C, sometimes it may be necessary to stir the water supply in the heat reservoir gently. As the apparatus functions, a temperature gradient is established in the heat reservoir, wherein the temperature at the surface is higher than at the bottom. To improve the output of the apparatus and lower the minimum temperature of the refrigeration chamber even further, raise the apparatus slightly and gently shake it back and forth (Photo 17). The agitation will homogenize the temperature in the heat reservoir and the temperature in the refrigeration chamber will diminish to attain a new minimum.


Photo 17

3.3.4 Studying fusion

It is also possible to study the fusion of the solid formed. The low thermal inertia of the refrigeration chamber and the small size of the sample being studied make it possible to observe the fusion of the solid formed in a few minutes after the power supply to the mini-freezer is turned off.

To do this, allow the system to warm up to the ambient temperature and measure the temperature as a function of time. Better results will be obtained if agitation is resumed as soon as there is a sufficient volume of liquid present in the refrigeration chamber to allow the user to work the agitator.

3.3.5 Working at elevated ambient temperatures

Under certain special environmental conditions (regions with hot climates, for example), the temperature of the local tap water may be too warm for the apparatus to function at its full output capacity. The chilling speed of the refrigeration unit and the minimum temperature attained by the apparatus may not be sufficient for obtaining satisfactory results.

If this is the case, the heat reservoir of the mini-freezer can be filled with water that has been prechilled in the refrigerator. A few hours before use, place a container of tap water in a refrigerator (1.5 L of water for each apparatus used). The chilling power of the apparatus can thus be increased considerably.

3.3.6 Very low temperatures

This apparatus can be used to produce the very low temperatures necessary for certain specific applications. To do so, fill the heat reservoir with an ice + water mixture at 0°C. By doing so, the refrigeration chamber will attain temperatures lower than -30 °C.

3.4 Manual measurements

3.4.1 Method

Monitoring the temperature manually is accomplished with a thermometer and a chronometer. The pairs of values T , t are recorded by the experimenter and the graph is drawn by hand on a sheet of graph paper (*Photo 18*).




Photo 18

3.4.2 Examples of results obtained

The graphs 19 and 20 shown below were created with the mini-freezer and a digital "Checktemp" thermometer. The temperature of the heat reservoir was 15°C at the start of the experiment.

Graph 19
 Study of the solidification of
 demineralized water
 Power supply: 12V

Graph 20
 Study of the solidification of
 cyclohexane
 Power supply: 12V

3.5 Computer-assisted measurements

The computer-assisted study of solidification can be performed with various types of materials.

3.5.1 Use of the VTT data acquisition console

Computer-assisted measurements are performed with a VTT data acquisition console equipped with a rapid response temperature sensor. The console can be used either in the autonomous mode or connected to a computer that uses the Compliss data acquisition software (*Photo 21*)

Photo 21

3.5.2 Use with the ESAO[®] Primo data acquisition console

Computer-assisted measurements are performed with an ESAO[®] Primo data acquisition console equipped with a rapid response temperature sensor. The ESAO[®] Primo data acquisition console is connected to a computer that uses the Compliss data acquisition software (*Photo 22*).

Photo 22

4 Storage and maintenance

4.1 Storage

4.1.1 Emptying the heat reservoir

Be sure to drain the heat reservoir before placing the mini-freezer in storage.

After you have disconnected the mini-freezer from the power source, open the fill hole (*Figure 1-8*) of the heat reservoir (as indicated in paragraph 2.2.1). Grasp the mini-freezer by the handle and tip it forward so that the water drains out through the fill hole.

4.1.2 Cleaning

The mini-freezer can be cleaned with a damp sponge and if necessary, household dishwashing detergent.

Be careful not to let water penetrate into the interior of the refrigeration device (if this does happen, refer to the instructions in paragraph 4.2.5).

Never use organic solvents (such as oxygenated, aromatic, or halogenated solvents), corrosive and oxidizing substances (sodium hypochlorite, etc.) to clean the mini-freezer. Doing so may result in permanent damage.

For long-term storage, clean the refrigeration chamber and the radiator of the mini-freezer and dry them carefully with absorbant tissue.

4.1.3 Storage

The mini-freezers can be stacked to facilitate storage and save cabinet space (*Figures 23 and 24*).

To stack, proceed as follows:

- Place the agitator in the refrigeration chamber (Paragraph 2.3).
- Return the thermometer holder to its storage position (Paragraph 2.2.2).
- Open the fill hole of the heat reservoir (Paragraph 2.2.1).
- Stack the mini-freezers on top of one another.


Figure 23


Figure 24

4.2 Care and maintenance

4.2.1 General information

This apparatus does not require periodic maintenance. However, it is recommended to follow the instructions below in order to maintain the long-term performance of the mini-freezer.

4.2.2 Cleaning the refrigeration chamber

To prevent dirt build-up and rusting, it is recommended to clean the refrigeration chamber carefully after each use. To do so, rinse the interior of the refrigeration chamber with demineralized water by filling it with the help of a pipette or a dropper. Then aspirate the rinse water and dry the inside of the refrigeration chamber with an absorbant tissue.

If necessary, the outside of the refrigeration chamber can be cleaned with a wet sponge and then dried with an absorbant tissue.

In the event of corrosion or the build-up of dirt deposits, gently scrub the inside surface of the refrigeration chamber with a piece of scouring pad for dishes until its metallic luster is restored.

4.2.3 Cleaning the heat reservoir tub

Calcium deposits may build up in the heat reservoir tub after several uses. It can be removed by pouring commercial vinegar or an aqueous solution of 5% acetic acid into the tub.

Proceed as follows:

- Separate the body of the apparatus and the heat reservoir tub (*figure 25*).
- Pour the vinegar or the dilute acetic acid solution into the tub of the heat reservoir.
- Wait until the effervescence stops.
- Rinse the tub with tap water and dry it with a piece of absorbant tissue.


Figure 25

4.2.4 Cleaning the radiator

As in the tub, calcium deposits may also build up on the radiator, mainly if the radiator remains immersed in water too long.

If the deposit is substantial, it may impede the efficacy of the heat exchanges between the radiator and the water in the heat reservoir. The radiator must then be cleaned.

Proceed as follows:

- Separate the cover and the heat reservoir tub.
- **Perform the following steps without replacing the cover so that the acid solution or the rinse water does not penetrate into the refrigeration device (figure 26).**
- Immerse the radiator in a container filled with commercial vinegar or an aqueous solution of 5% acetic acid.
- Allow it to react for two to three minutes.
- Rinse the radiator copiously by immersing it in a container filled with tap water.
- If necessary, repeat these last three steps until the calcium deposit has been eliminated.
- Dry the radiator with a piece of absorbant tissue.

Never clean the radiator with a concentrated acid solution. Doing so could damage it permanently!


Schéma 26

1 – Body of the refrigeration device.
2 – Radiator.

3 – Container (wash bowl or cooler).
4 – Commercial vinegar or 5% acetic acid solution

4.2.5 Drying the refrigeration device

This operation should only be performed if a significant amount of water has penetrated into the refrigeration device.


Figure 27

Remove the protective cover from the refrigeration device as indicated in figure 27.

Wipe up the liquid by gently blotting the affected area with a piece of absorbant tissue.

Be very careful not to damage the wiring of the device.

Replace the protective cover on the refrigeration device.

5 Accessories

We recommend the following accessories in order to exploit the performances of the mini-freezer to their fullest potential.

5.1.1 Power supply

Name	Item no.
Evolution® F6F12 / 5 A transformer	281 083
Red/black safety cables, 50 cm in length	283 357

5.1.2 Manual measurements

Name	Item no.
Checktemp -50/-150°C digital thermometer	251 030
1/100 ^e s digital chronometer	351 037

5.1.3 Computer-assisted measurements

VTT® console option

Name	Item no.
VTT console	471 000
Rapid response thermometric sensor (-20/+120°C)	472 053
Compliss VI software	000 512

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Name	Item no.
ESAO® Primo console	451 038
Rapid response thermometric sensor (-20/+120°C)	472 053
Compliss VI software	000 512

6 Troubleshooting

6.1 Correct function indicator

The mini-freezer is equipped with an indicator light (*Figure 1-2*). When it is on, it shows the user that the device is functioning normally.

When the indicator light is off, this means that the thermal circuit breaker has been tripped or even that the apparatus has malfunctioned. In the event that this occurs, refer to paragraph 6.3.

6.2 Safety features of the apparatus

6.2.1 Polarity

The refrigeration device is protected against polarity inversions. If the power supply to the apparatus is not hooked up in accordance with the indicated polarity, the circuit is cut off and the function indicator light does not come on.

6.2.2 Overheating protection

The refrigeration device is protected against all abnormal temperature increases. A thermal circuit breaker will cut off the power supply in the event that the device overheats.

The usual reasons why the thermal circuit breaker trips are:

- Continued operation of the apparatus after the minimum temperature in the refrigeration chamber has been attained. This is most likely to occur at high ambient temperatures.
- The lack of water or too low of a water level in the heat reservoir.
- The temperature of the water in the heat reservoir is too high.

6.3 Finding the cause of and fixing malfunctions

6.3.1 The function indicator light does not come on

Check the apparatus over in the following order:

- Check the power supply and make sure that the current is on.
- Check the connection between the mini-freezer and the power supply, as well as the polarity of the connection.

If you are not able to identify the cause of the malfunction at this stage, it is possible that the problem has caused the thermal circuit breaker to trip. In this case, proceed as follows:

- Disconnect the mini-freezer from its power supply.
- Make sure that there is sufficient water in the thermal reservoir.
- If it is empty, fill it up with **cool water** (Refer to the instructions in paragraph 2.2.1).
- If water is already present in the heat reservoir, empty it and refill it with **cool water**.
- Allow the apparatus to cool for a few minutes.
- Reconnect the mini-freezer to its power supply, and the function indicator light should come back on.
- If the thermal circuit breaker trips again soon afterwards, unplug the apparatus and allow it to cool for a few more minutes.

If the apparatus is still not working after performing these operations, it is because the refrigeration device itself is malfunctioning. **Contact our Technical Support.**

7 After-sales service

The apparatus is under a 2 year warranty. In case of defects, the material should be returned to our shops.

For all repairs, adjustments, or replacement parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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Service gratuit * :
+ 33 (0)2 32 29 40 50

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Rubrique FAQ

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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