Appareil de chauffage

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Heater

Laboratory electric spout LAB2

Réf :
701 280

Français – p 1

English – p 8

Version : 8107

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Burner electric LAB 2

1. Description

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
1. Description

The LAB 2 laboratory electric spout is a heater that can be substituted for Bunsen burner to heat test tubes, containers and to conduct several chemistry experiments (Iron and sulfur reaction, charcoal inflammation...).

- 1 – Support plate for flat bottom containers (Beakers, balloons...)
- 2 – Thermal protection cowling
- 3 – Heating the resistance indicator
- 4 – Timed adjustment from 10 to 650 °C
- 5 – Protection fuse
- 6 – Power On indicator

2. Specifications

Weight : 3.30 kg
Dimensions (mm) : L. 170 x l. 130 x h. 195


- Power: 500 W
- Power supply: 230 V / 50 Hz (2P+T)
- Rapid fuse protection: F. 5 A.H.250 V
- Power On indicator: Red
- Temperature indicator: Blue
- Temperature adjustable from 10 to 650°C
- Integrated controller
- Protection from overheat

Environmental conditions:

- Maximum use room temperature: 40°C
- Relative humidity below 80 %

Pull the power plug which has to remain accessible to power down the device.

The laboratory electric spout LAB 2 complies with the standards EN 61010-1 and EN 61010-2-010.

3. Use

3.1 Switching On and Off

- Place the device on a stable and plane surface
- The device will be connected to a network corresponding to its electrical specifications with earth connection and differential protection.
- After connection to the network, the power On indicator turns on (red indicator).
- Rotate the adjustment knob from 0 to MAX, depending on the required temperature. The blue temperature indicator remains lit until obtaining the desired temperature, and then it turns off (except for the MAX position).
- To turn off the device: reset the adjustment knob to position 0. The red indicator turns off. The device is completely off.

3.2 Heating Containers

To heat flat bottom containers (balloons, beakers, etc...): place them directly on the support plate.

To set the round bottomed containers on the LAB 2 electric spout, or to ensure their secure fastening, we suggest using the support for glassware (ref. 703 522) or the rod for electric spouts LAB 2 and LAB 3 (ref. 703 087).

Round-bottomed flask retention with the support for glassware (ref. 703 522)

Setting a round-bottomed flask using a rod for electrical spouts LAB 2 and LAB 3 (ref. 703 087).

Warning

The round bottom (of flasks, beakers, etc...) shall not, in any case, be higher in surface than the protection cowl, being a maximum of 11 cm.

3.3 Heating Test Tubes

To heat the test tubes, manually hold the tube using a wooden clip or place it in the glassware support (ref. 703 522), in contact with the grate.

Manually heating a test tube.

Setting a test tube using a support for glassware (ref. 703 522)

For heating without **flashing out** the liquid in a test tube, set the adjustment knob to position MAX so that the air column wraps the full tube length.
If the device is cold, before each use, set the temperature to MAX for 45 seconds in order to create the heat column, and then set it to the desired temperature.

Warning


In order to avoid flashing out, place the tube at the center of the device on the warm air column without shaking it.

3.4 Producing a Flame

This device allows conducting chemistry experiments that require using a flame (combustion of metallic powder...).

For this, use a spatula to put an Ecoflam tab (ref. 101 012) on the LAB 2 electric spout grid, the resistance of the latter having been previously preheated. The drop ignites spontaneously and produces a flame that can be used for about 1 minute.

The combustion reactions can be conducted directly above the electric spout grid.

Production of a flame on the electric spout grid with an Ecoflam tab

Bottle of Ecoflam tabs (ref. 101 012) containing about 60 tabs

3.5 Combustions

The heat flux produced by this device allows igniting several materials such as steel wool, a charcoal stick or a magnesium ribbon. Once the device is preheated, hold the material to be ignited using a wood clip or place it above the flow of hot air (centered on the grid). Wait for a few seconds until the combustion reaction starts.

3.6 Handlings Plates

They are used to conduct combustion reactions with a flame or pyrolysis, and dehydration reactions or reactions requiring a high temperature (reaction between iron and sulfur). These plates allow accommodating combustion or reactions products.

Using the combustions plate for reactions requiring a flame.



Plate for reaction and combustion of the mixture Iron/Sulfur (Ref. 703 010).


Plate for combustion of metallic powders (Ref. 703 015).

3.7 Handlings Plates for Microbiology

These accessories allow using the LAB 2 electric spout as part of microbiology handlings in a sterile area.

In fact, these plates allow operating the sterility cone generated by the device. Just place the culture media on the latter once fixed to the casing of the device in order for the handlings to be carried out in sterile area.

Handling plate for microbiology (Ref.: 701 299).

4. Safety

4.1 Risk of Burns

- The LAB 2 laboratory electric spout is a heater that should be handled carefully.
- It is not advisable to move the device when it is heated or turned on.

4.2 Personal Safety

- Do not carry out any experiment before being previously aware of the risks related to the chemical products used.
- Work with the personal and collective safety equipment adapted to the risks specific to the chemical products used (Hood, Safety goggles...).
- do not use this device to heat flammable liquids, products or mixtures of products with risks of inflammation or explosive decomposition.

4.3 Gas Emissions

Work in a ventilated room if there is a risk of gas emissions during the handlings carried out.

4.4 Use Precautions

Do not leave the device unattended in continuous operation.

4.5 Work Surface

- Never use the laboratory electric spout on an inflammable surface.
- Never work at less than 30 cm from a wall.

4.6 Warning

JEULIN Company alerts the user that the safety of this device could be undermined if used in a manner not specified in this manual.

5. Care and Maintenance

- This device does not require periodic or preventative maintenance. In case of pollution, it should be cleaned by taking the following precautions:
 - **the device should be disconnected from the power supply.**
 - cleaning (regardless of the product used) should only be made if the device is cold.
- To clean the inside of the heating element from the combustion residues that could be in it, simply turn around the device unplugged in advance and the cold resistance and shake vigorously. Residues will fall away on their own.
- This device does not require maintenance of the electrical part.

- **In case of failure**

If the power up red indicator does not turn on:

- Set the adjustment knob to 0.
- **Unplug the mains plug.**
- Completely turn the device around.
- Change the fuse.

5 ampere fast-blow fuse. High breaking capacity 250 volts (compatible CEI 127).

If, after verifying the fuse, the device still does not operate, contact our company:

This device must in no way be dismantled by the used for any reason. Any intervention should be carried out by JEULIN Company or its representative.

6. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux