Appareils de laboratoire

Chauffage

Laboratory equipment

Heating

Ref :
701 323

Français – p 1

English – p 9

Version : 7001

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LAB 3 Laboratory electric burner

1 Description

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
NOTES

1 Description

The laboratory electric burner LAB 3 is a heating apparatus that can be used as a substitute for a Bunsen burner for heating test tubes and containers, and for performing diverse chemistry experiments (iron and sulphur reactions, charcoal combustion...)

1 - Support grill for flat-bottomed containers (beakers, flasks...)
2 - Protective housing that always remains cool
3 - Temperature indicator


4 - Time delayed temperature control with a range of 10 to 750 °C
5 - Safety fuse
6 - On/off indicator light

2 Specifications

Weight: 3.45 kg

Dimensions (mm) **l.** 130 x **w.** 170 x **h.** 195

- Power output: 900 W
- Voltage/frequency: 230 V - 50 Hz (2 P+T)
- Fast-acting safety fuse: F. 6.3 A.H.250 V
- On/off indicator light: red
- Temperature indicator light: blue
- Temperature adjustable from 70 to 750 °C.

- Built-in thermostat
- Protection against overheating

Temperatures corresponding to the LAB 3 graduations:

Graduation	Temperature (°C)
1	70
2	210
3	280
4	420
5	490
6	630
7	750

Environmental conditions:

- Maximum permissible ambient temperature: 40 °C
- Permissible relative humidity: less than 80 %

The power to the apparatus is turned off by unplugging it; therefore, the plug should always be accessible.

The laboratory electric burner LAB 3 complies with the EN 61010-1 and EN 61010-2-010 standards.

3 Use

3.1 Power on / off

- Place the apparatus on a stable, level surface.
- The apparatus must be plugged into a circuit that meets its electric requirements, with grounding and differential protection.
- After plugging in, the on/off indicator light comes on (the red light)
- Turn the control knob from 0 to MAX, according to the temperature required. The blue temperature indicator light will stay on until the desired temperature is reached, then it goes out (except when the knob is set on MAX).
- To turn the apparatus off: set the control knob back to 0. The red indicator light will go out. The apparatus is now turned off completely.

3.2 Heating containers

- To heat flat-bottomed containers (beakers, flasks...): place them directly on the support grill.
- To place round-bottomed flasks on the LAB 3 electric burner or to enable a safer maintaining of flat-bottom containers, it is recommended to use the glass holder (Item no. 703 522), or the grid for LAB 2 and LAB 3 electric burners (réf. 703 087).

Maintaining the round-bottomed flask using the glass support (ref. 703 522)

Fixation of a round bottomed-flask using a rod for LAB 2 and LAB3 electric burners. (ref. 703 087).

CAUTION:

The surface area of the flat bottomed containers (flasks, beakers, etc..) must not overlap that of the protective housing, i.e., they must be no greater than 11 cm wide.

3.3 Heating test tubes

For test tubes: hold the tube with a wooden clamp or place it in the glass holder (Item no. 703 522), in contact with the support grill.

Manual heating of a test tube

Fixation of a test tube using a glass support (ref. 703 522)

In order to **avoid boiling over** when heating liquids in a test tube, set the control knob to the MAX position so that the air column envelops the entire length of the tube

If the apparatus is cold: before each use, set the temperature to MAX for 45 seconds to create a column of heat, then set it to the desired temperature.

IMPORTANT:

To avoid boiling over, either place the tube in its holder or avoid shaking it.

3.4 Producing a flame

This device enables you to perform experiments in chemistry requiring the use of a flame (metal powder combustion ...).

In this aim, using a spatula, place an Ecoflam pellet (ref. 101 012) on the grid of LAB 3 burner which resistance has been first pre-heated. The pellet spontaneously ignites and produces a flame which can be used during about 1 mn.

The reactions of combustion can be directly performed above the grid of electric burner.

Production of a flame on the electric burner with an Ecoflam pellet.

Box of Ecoflam pellets (ref. 101 012) containing about 60 pellet


3.5 Combustions

The heat flow produced by this device enables the ignition of diverse materials such as steel wool, a rod of charcoal or a magnesium band. Once the device has been preheated, take the material to be ignited using a wood clamp and place it above the warm air flow (centered on the grid). Wait a few second until the reaction of combustion begins.

3.6 Plates for experiments

To perform reactions of combustion with flame or pyrolysis reactions, dehydration or reactions requiring a high temperature (reaction between iron and sulphur). These plates enable the collection of combustion products and reactions products.

Using of the combustion plate during reactions requiring a flame. ➔

↑
 Plate for reaction and combustion of iron / sulphur mixing (Ref. 703 010).

↑
 Plate for the combustion of metal powder (Ref. 703 015).

3.7 Plates for experiments in microbiology

These accessories enables the use of LAB 3 electric burner in the frame of microbiology experiments in sterile area.

Indeed, these plates enable the exploitation of sterility cone produced by the device. You only have to place culture medium on them once they are fixed on the device's frame so that the experiments are performed in a sterile area.

Plates for experiments in microbiology (Ref. 701 299).

4 Safety

4.1 Burn danger

- The laboratory electric burner LAB 3 must be handled with caution.
- It is not advisable to move the apparatus while it is heating or turned on.

4.2 Personal protection

- Do not perform any experiment without first knowing the risks associated with the chemicals used.
- Use the proper personal protective and laboratory safety equipment (fume hoods, protective goggles...) for the type of chemicals being used.
- Do not use this apparatus to heat inflammable liquids, chemicals or mixtures of chemicals that pose a fire or explosion hazard.

4.3 Fumes

Work in a ventilated room if there is a danger of fumes being generated from a given procedure.

4.4 Use precautions

Never leave the apparatus running unattended.

4.5 Work surface

- Never use the apparatus on an inflammable surface.
- Never work closer than 30 cm to a wall.

4.6 Warning

The JEULIN Company warns the operator that the safety of this apparatus may be compromised if it is used in any manner not specified in these instructions.

5 Care and maintenance

- This apparatus needs no periodic or preventative maintenance. In the event of contamination, this apparatus must be cleaned while observing the following precautions:
 - **the apparatus must be unplugged.**
 - do not clean (no matter what product is used) until the apparatus has cooled completely.
- The electrical components of this apparatus require no maintenance.
- **If the apparatus is not working:**
If the red indicator light does not come on:
 - Set the control knob to 0.
 - **Unplug the apparatus.**
 - Turn the entire apparatus over.
 - Change the fuse.

Fast-acting 6.3 ampere fuse. Maximum capacity 250 volts (in compliance with CEI 127).

If the apparatus still does not work after checking the fuse, contact us.

Under no circumstances should the operator dismantle this apparatus. Such operations should be only performed by the JEULIN Company or one of their representatives.

6 After-Sales Service

This material is under a three year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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Service gratuit * :
+ 33 (0)2 32 29 40 50

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Rubrique FAQ

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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