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# Chemical Measurement

Colorimetry

Réf :  
701 648

Français – p 1

English – p 5

Version : 8107

**Colorimètre Color 1G-7**  
***Color 1G-7 colorimeter***

## 1. Description

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**0 825 563 563\***

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## 1. Description

The Color 1G7 colorimeter can measure the light absorbance and transmission through a colored solution over 7 wavelengths. It may be used to follow enzymatic or chemical kinetics and to perform a colorimetric assay (Beer-Lambert's law).


Figure 3: View from above



Figure 4: Rear view

## 1.1 Spectral Range

The Color 1G7 colorimeter is equipped with 7 light-emitting diodes from 470 nm to 660 nm.

470      525      570      590      605      626      660

## 1.2 Navigation

The central wheel is used for the navigation. Switching from one tab to another is done by rotation and validation with a press.

## 2. How to setup

- Connect the Color 1G7 colorimeter through the mini USB plug.
- Place a tank containing the solvent used in the tank location, and firmly press the tank.

Note: make sure you respect the direction of the tank. The grooves must be placed on both sides as shown on the Color 1G7 colorimeter diagram

- Cover the tank.
- Use the wheel to select the wavelength. To do so, place the dotted rectangle above the required wavelength and validate by pressing on the center of the wheel. Once selected, the rectangle turns black.
- Select the measure type: absorbance or transmission. For this, surround "Absorbance" or "Transmission" with a dotted rectangle using the wheel and validate your choice by pressing the center of the wheel. Your selection turns black.
- Briefly press the selected measure type to proceed with the calibration and press the center of the wheel to validate the calibration.
- Then place the tank containing the solution to be analyzed. The result (absorbance or transmission) directly appears on the screen.

### 3. The Zero Adjustment

The zero adjustment is necessary before each measurement or series of measurements over a wavelength. This step allows defining the reference light intensity, i.e. the light intensity received in the presence of the solvent.

- Place a tank filled with solvent into the tank carrier
- Use the wheel to place the dotted rectangle on the selected measure type (absorbance or transmission)
- Press on the center of the wheel
- Select **Calibration**
- The displayed value thus indicates 100% for the transmission or 0 for the absorbance.

This operation is necessary for each change of wavelength. Moreover, a 5 minutes warm-up time is recommended before each zero adjustment.

### 4. Analog Output

The Color 1G7 colorimeter is equipped with two banana sockets for digitalizing measurements through a suitable data logger. By way of example, we will use the FOXY® interface.

- Use a banana cord to connect the Color 1G7 yellow terminal to the red terminal of a direct input on the Foxy
- Repeat this action with the black terminal
- In the software, place the direct input icon on Y axis
- Select  $\pm 5V$  in the Caliber tab
- In the Customize tab, associate a name and a unit to the new value.
- According to the chosen measure (absorbance or transmission), fill in the correlation table as follows:

Caliber	Quantity	Measure	Color	Customize
	Quantity Abs			
	Unit Abs			
	Point n°1 0	0		
	Point n°2 3	3		

*Absorbance*

Caliber	Quantity	Measure	Color	Customize
	Quantity T			
	Unit %			
	Point n°1 0	0		
	Point n°2 100	1		

*Transmission*

Note: These values are listed on the back of the Color 1G7 device.

- Click on the floppy disk to save the settings 

A zero setting on the Color 1G7 device is required when changing the LED (see previous point). The direct inputs adjustment does not change.

## 5. Technical Specifications

<b>Power supply</b>	5V, 1A via mini-USB (supplied)
<b>Sizes</b>	160 x120 x70 mm
<b>Light source</b>	Light emitting diode
<b>Wavelength</b>	470, 525, 570, 590, 605, 626, 660 nm
<b>Absorbance</b>	-3 to 3 abs
<b>Transmission</b>	0 to 100 %
<b>Analog outputs</b>	Ø 4 mm safety banana sockets -3V/3V corresponding to -3 abs/3abs of absorbance 0V/1V corresponding to 0/100 % of transmission

## 6. After-sale Services

2 years warranty.

For any repair, adjustment or spare parts concerning this device during or after the warranty, kindly contact:

**JEULIN – TECHNICAL SUPPORT**  
**468 rue Jacques Monod**  
**CS 21900**  
**27019 EVREUX CEDEX France**

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

### Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux