Mesure

Initio 2

Measure

Initio 2

Réf :
701 653

Français – p 1

English – p 9

Version : 7105

Conductimètre Initio 2
Initio 2 conductivity meter

SOMMAIRE

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
1. Description

The Initio 2 conductivity meter is an instrument that measures the ionic conductivity of a solution making use of a conductivity measurement cell. It can also save and store data, and restore them later onto a computer.

A WIFI module (ref 451 470) is also available to display and save the measurement on a tablet or a PC with WIFI.

The Initio 2 conductivity meter is fitted with an analogue output and can be used with an ExAO interface with direct inputs through banana connectors or a Groomy® interface using the jack.

2. Use

It is recommended that the Initio 2 battery should be fully recharged when it is used for the first time. This should be done using the μ USB charger (7) delivered in a separate package. In general, it is recommended that complete charge/discharge cycles should be carried out (wait until the instrument is fully discharged before completely recharging it).

2.1 Start up

- Press the power on button (1) on the left side of the instrument.

The measurement is displayed directly on the instrument screen.

- To switch the Initio 2 off, do a short press on the power on button (1) and validate pressing the knob (3).

2.2 Browsing

Initio 2 has three buttons for browsing :

Selection-validation (3):

The central knob (3) is used to browse through the Initio 2 menus. Do a brief press on the OK button on the middle of the knob to validate.

Return (2):

The "Return" button (2) changes from one type of display to the other from the main screen in the following order:

- Measurement
- Table
- Graph

Within a particular menu, it acts as a Back button.

Menu (4):

This button displays the Initio 2 main menu.

2.3 Range

The Initio 2 Conductivity meter has two ranges

- 2 mS/cm
- 20 mS/cm

To change from one to another:

- Go to the main screen (measurement display) using the Return button
- Press the validation knob
- Select the required range
- Validate by pressing the knob

2.4 Use with the CTM probe

The CTM probe does not have a built-in temperature sensor. However, the measurement can be offset making use of an external thermometer when the calibration is made:

With external thermometer

- Press the Menu "button"
- Select Calibration
- Enter the temperature read with the thermometer using the rotary knob
- Select Calibrate
- Enter the value of the standard solution at 25°C and validate

Without external thermometer

- Press the Menu "button"
- Select Calibrate
- Select Without offsetting

- Select Calibrate
- Enter the standard value and validate

2.5 Use with the CTA probe

The CTA probe has a temperature sensor that automatically offsets the measurement made by the cell. The calibration is made with the value of the standard solution at 25°C. However, the offsetting can be deactivated (display the effect of temperature on the conductivity of the solution)

- Press the Menu button
- Select Calibrate
- Select With or Without offsetting
- Select Calibrate
- Enter the value of the standard solution at 25°C and validate

2.6 Data acquisition

The Initio 2 has an internal memory to save measurements. There are two available modes:

- Manual acquisition

Acquisition is made point by point. The user saves each measurement manually (see below).

- Time acquisition

Acquisition is made as a function of time, for a duration input by the user.

2.6.1 Manual acquisition

- Press the Menu button (4).
- Use the knob (3), to choose **Acquisition** and validate by pressing briefly on the knob (3).
- Select **Manual**
- To save a measurement, choose the **Save** icon (white circle) and press knob (3) briefly to validate.
- When you want to stop acquisition, select the **Stop** icon (black square) and validate by pressing knob (3). Validate the choice in the dialogue box.

2.6.2 Time acquisition

- Press the Menu "button" (4).
- Use the knob (3) to choose **Acquisition** and press the knob (3) briefly to validate
- Select **Time**.
- Use the knob (3) to select the number defining the duration and press the knob (3).
- Modify this duration by turning the knob (3) and then validate by pressing the centre of the knob.
- For a loop acquisition (restart the acquisition at zero after the defined time) select **In loop**.
- Validate the final choice by selecting **OK**.

Acquisitions as a function of time last between 60 seconds and 50 hours. The interval between two measurements is adjusted automatically.

2.6.3 Save data

Once the acquisition is complete (time or manual), you need to save it in order to retrieve points on a computer.

- Press the Menu "button" (4)
- Select **File** and validate
- Select **Save**

You can only save one acquisition. Select **Open** to open the most recently saved acquisition.

2.6.4 Transfer data to a computer

Once the acquisition is complete, you can transfer the data onto a computer for processing:

- Install the Initio 2 software available at address www.jeulin.com/Initio 2
- Connect the Initio 2 to the computer using the USB- μ USB cable
- Open the Initio 2 software
- Load the saved acquisition by clicking on the icon:

Data are then displayed on the screen. Do a right click on the graph to access the different graphic tools.

3. Data display

You can display data in the form of a table or a graph at any time. Simply press the Return key (2). You thus change in sequence from the measurement to the table and then to the graph. You can continue the acquisition regardless of what is displayed.

The table or graph will be displayed in full screen after 5 seconds without taking any action. Turn the knob if you want to obtain the associated functions once again.

The arrows at the bottom right in the table are used to move through the rows of the table.

4. Analogue outputs


The Initio 2 has analogue outputs on banana connectors and on the 3.5 mm jack.

4.1 Banana connectors


- Connect the Initio 2 to an ExAO interface provided with direct inputs using banana cables
- Configure the direct input depending on the range:
 - o 0-2ms/cm
 - 0 mS/cm corresponds to -5V
 - 2 mS/cm corresponds to +5V
 - o
 - 0 mS/cm corresponds to -5V
 - 20 mS/cm corresponds to +5V

Example with the Foxy® console

- Start the Generalist
- Put the icon for a direct input  on the required axis
- Select the $\pm 5V$ range

- Fill in the customisation tab as a function of the range, for example:

- Click on the diskette to save the parameter settings 

The parameter settings are complete, the displayed magnitude is now the conductivity with symbol σ and unit is mS/cm.

4.2 Jack

The 3.5mm jack is used for example to connect the Initio 2 to a Groomy® interface (ref 181 000). The output voltage varies between 0 and 5V. We have:

- -20°C corresponds to 0V
- 120°C corresponds to 5V

5. WIFI communication

The WIFI module (ref 451 470) is used to display Initio 2 data, for example on a tablet.

5.1 Installing the WIFI module

The WIFI module is installed when the Initio 2 is switched off.

- Raise the Initio 2 stand
- Remove the back by sliding it downwards
- Put the WIFI module into place, replacing the back
- Push the module all the way into place until it click fits. It must be possible to retract the stand without forcing

5.2 Display on tablet

Once the WIFI module has been installed, the tablet has to be configured:

- Press the Menu "button"
- Select Parameters

- Select Network

Information about the module is displayed on the screen.

- On the tablet, choose the name of the WIFI network displayed on the Initio 2.
- Open an internet browser
- Enter the address displayed by the Initio 2 in the browser address bar.

6. Parameters

The various information on the Initio is accessed by pressing Menu and selecting Parameter → Information

This page displays the serial number of the instrument and the software version, the level of the battery charge and the working language. The working language can be modified using the knob.

7. Update

On our website at the address <http://en.jeulin.fr/en-resources/english-resources/software/free-software.html>, the last update for your Initio® 2 is available. Its version is indicated in the designation.

In order to know the version of your Initio® 2, press the “Menu” button and then “Setup” “Informations”.

8. Technical characteristics

Display: Backlit graphic screen, resolution 128 x 64 (Visible size: 66.5 x 33.2mm)

Range: 2 mS.cm⁻¹ and 20 mS.cm⁻¹

Measurement range:

0 to 2 mS.cm⁻¹; resolution 0.001 mS.cm⁻¹ ; precision ±1% of the full scale
(Adapted standard solution : 1413 µS.cm⁻¹ ref. 106 191)

0 to 20 mS.cm⁻¹; resolution 0.01 mS.cm⁻¹ ; precision ±1% of the full scale
(Adapted standard solution : 12 880 µS.cm⁻¹ ref. 106 206)

Analogue output:

-5 V to +5 V on 4 mm safety banana sockets
0 to +5 V on 3.5 mm Jack

Calibration: At one point

Operating mode: CTA or CTM

- CTA mode: Compatible with the CTA probe (ref. 701 330) on mini DIN plug
- CTM mode: Compatible with the CTM probe (ref. 701 639) on banana connectors

CTA and CTM probes: To be ordered separately

Automatic temperature offsetting: 0 to 50°C

Optimum offsetting range: 15 to 35°C

Communication: USB, WIFI optional

Dimensions: 95 x 123 x 44 mm

9. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

+33 (0)2 32 29 40 50