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Vaporisation FR

# ***Laboratory apparatus GB***

*Vaporization GB*

Réf :  
701 679

Français – p 1

English – p 3

Version : 8105

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***Mini-vaporizer GB***

## 1. Description

### 1.1 Composition

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
### 1.1 Composition

The apparatus (réf 701679) consists of:

- 1 transparent low volume vial
- 1 heating resistor suitable for handlings starting from middle school
- 1 stable support to avoid spillage

### 1.2 Heating Resistor

The heating resistor is designed to be used by students starting from middle school: its upper part doesn't heat to avoid burns during the handlings.

### 1.3 Support

The support, highly stable, is fitted with 2 banana sockets located at the back for the heating resistor power.


This product must be operated with a power: 12V – 2A mini.  
*We recommend a power: réf 281483.*

### 1.4 Compatible Thermometers

The mini-vaporizer accepts all thermometers with a 5 mm diameter temperature sensor.

For the Jeulin Computer Assisted Experimentation thermometers, a rubber element controls perfectly the sensor positioning in the liquid volume.

For any other thermometer, remove the rubber element to give free passage to sensors of higher diameters.


## 2. Achievable Experiments

This product is dedicated to the study of the process of changing from a liquid to a gas. The experimenter will be able to test different liquids such as: water, salt water, alcohol, or others.

Warning: in order to compare the vaporization time, the same weight should be respected for the different liquids.

## 3. Result Example: Case of Water

Equipment used:

- Power supply: réf 281483
- Thermometer sensor: réf 482204
- Acquisition console: Foxy®, Tooxy® or AirNeXT®
- Mini-vaporizer: réf 701679

## 4. After-Sale Services

The device is under a 2-year guarantee, it must be sent back to our workshops. For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
468 rue Jacques Monod  
CS 21900  
27019 EVREUX CEDEX FRANCE

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

### Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux