**TP de chimie**

***Chemistry tutorial***

**Ref :  
702 081**

**Français – p 1**

**English – p 3**

**Version : 6009**

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This set is designed for the teacher. It can be used for representing the close-packed or split model of different molecules and for solving the problems in stereochemistry.

## 1 Composition

Qty.	Element		Colour	No. of holes / type	Ø mm
24 6 2 6	Carbon	C	Black	4 / tetrahedron $sp^3$ 3 / trigonal 2 / linear 5 / tripyramid	23
40	Hydrogen	H	White	1	17
1 12	Oxygen	O	Red	4 / tetrahedron 2 angular	23
1 4	Nitrogen	N	Blue	4 / tetrahedron 3 / tetrahedron	23
1 1	Sulphur	S	Yellow	4 / tetrahedron 2 / angular	23
12	Halogen		Green	1	23
4	Phosphorous	P	Purple	4 / tetrahedron	23
2 1	Metal		Grey	1 2 angular	17 23
3	Orbital		Purple	Length : 39 mm	
55	Medium link		Grey		19/31
25	Long flexible link		Grey		32/43
60	Short link		White		2/10
1	Dismantling tool				

## 2 Types of links

The grey links are used for symbolising the simple covalent links.

The long flexible grey links allow representing the multiple covalent links in the split models (double or triple links).

The short white links are used instead of the medium size links for creating close-packed models.

## 3 Notes on the carbon atoms

- The atom  $sp$  linear 2 holes  $180^\circ$  is used for building linear models (e.g. acetylene).
- The atom  $sp^2$  trigonal 3 holes  $120^\circ$  is used for the molecules placed on a plane (e.g. ethylene, benzene).

- The atom  $sp^3$  trigonal dipyramid 5 hole has two holes arranged as in the atom 2 hole linear  $sp$  and three holes arranged at  $120^\circ$  as in the atom  $sp^2$ . It is therefore an element that allows creating a large number of structures.

## 4 Using the dismantling tool


- Place the tool on the link to be removed (the face with the markings should be turned upwards).
- Push the key horizontally and completely under the link which will be lifted by 2 mm. Do not try to lift the tool before positioning it correctly; this will create unnecessary tensions on the edges of the dismantling tool.
- Release the link gently through a lifting movement on the tool while holding it with the thumb.

## 5 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**+33 (0)2 32 29 40 50**

*The Molymod<sup>®</sup> molecular models are made in England by SPIRING Enterprises Ltd.+*

*The registered trademarks MOLYMOD<sup>®</sup>, MOLYDOME<sup>™</sup> and ATOMOD<sup>™</sup> are the exclusive properties of Spiring Enterprises Ltd and all rights are reserved.*

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*Rubrique FAQ*

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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