Chimie structurale

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Structural Chemistry

Properties of Matter

Réf :
702 852

Français – p 1

English – p 6

Version : 8108

Structure de l'atome

Atom Structure

1. Description

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1. Description

1.1 Presentation

This set of tokens is used for modelling the atom and its internal structure: atomic nucleus (nucleons: protons, neutrons), electrons.

Through easy visual handlings, it is possible to:

- distinguish atom and element
- understand the composition of the atomic nucleus by distinguishing protons and neutrons;
- understand the electrical neutrality of an atom;
- understand the monoatomic ions formation;
- understand the net structure of atoms, thus of the matter.

1.2 Composition

The box contains


- 20 red tokens:
- 20 blue tokens:
- 20 white tokens:

Protons
Neutrons
Electrons

The number of tokens in the box allows modeling the atoms up to $Z = 20$ (element Calcium).

- ⇒ The label has a sidebar to be filled with the exact weight of the whole box.
This tip allows to quickly check if it is full or not at the end of the PW session.


The board for atomic model is available under A4 format to be downloaded alone or found at the end of this manual:

2. Experiments


2.1 Production of Atomic Models

2.1.1 The Simplest Atom of the Universe – Hydrogen H


Hydrogen model 1: ^1H
Composition: 1 proton + 1 electron

2.1.2 Two Different Atoms of the Same Element – Helium He

Helium atom 4: ^4He
Helium element: **He**


Helium atom 3: ^3He
Helium element: **He**

Both atoms have the same number of protons and electrons ($Z = 2$) but the helium ^4He has 2 neutrons while the ^3He only has one.

2.2 Production of Monoatomic Ionic Models

2.2.1 Positive Ion = Sodium Example


The nucleus common to both particles has 11 protons (red tokens) and 12 neutrons (blue tokens).

The sodium atom being neutral, it has as many electrons as protons, which means 11 electrons.

The sodium ion, which is positive, has 1 less electron, meaning 10 electrons.

2.2.2 Negative Ion = Chlorine Example


The nucleus common to both particles has 17 protons (red tokens) and 18 neutrons (blue tokens).

The chlorine atom being neutral, it has as many electrons as protons, which means 17 electrons.

The chlorine ion, which is negative, has 1 extra electron, meaning 18 electrons.

3. Board to Print

NUCLEUS

ATOM

4. After-Sale Services

The warranty is of 2 years.

For all settings, contact out **Technical Support** on **0 825 563 563**.

The material must be returned to our workshops and for all repairs or spare parts, kindly contact:

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux