Chimie

Chemistry

**Ref :
723 032**

Français – p 1

English – p 5

Version : 6010

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NOTES

1 Description

This chromatography vessel is especially designed for an introduction to paper chromatography. Constructed of unbreakable polypropylene, it is transparent in order to allow the student to easily observe the phenomena taking place during the elution. The bottom of the vessel is equipped with a mark to indicate the level of the eluant for easy filling. The lid is equipped with a device to secure the chromatography paper.

Photo 1

- 1 – Unbreakable polypropylene vessel.
- 2 – Eluant level mark.
- 3 – Polyethylene lid.
- 4 – System for securing the chromatography paper.

2 Use precautions

This vessel is designed for an introduction to paper chromatography using an aqueous eluant. Aqueous saline solutions (water + salt) and aqueous mixtures (water/ethanol, water ethanoic acid, etc.) may also be used.

Under no circumstances is this vessel to be used with organic solvent-based eluants (halogens, aromatics, alkanes, esters, etc.). Doing so may damage the vessel and the safety of the user is not guaranteed.

3 Use

3.1 Preparing the vessel

1. Remove the lid from the vessel.
2. Fill the vessel up to the level mark with the appropriate eluant.
3. Replace the lid on the vessel.

3.2 Inserting the paper

1. Open the chromatography paper support system by pressing on the lid on either side of the slot.
2. Insert the sheet of chromatography paper.
3. Slide the sheet slowly until it is dipped in the eluant (Figure 2). Be careful not to move the vessel when doing so to avoid making ripples in the eluant.
4. Release the pressure on the lid so that it closes and pinches the chromatography paper.


Schéma 2

3.3 Tricks and hints

This vessel is designed for paper chromatography. It is recommended to use chromatography paper; however, ordinary filter paper may work for certain applications.

The unique property of chromatography paper is that it is made of cellulose fibers aligned in a preferred direction. It thus has a strong “capillary force” that favours the upward movement of the eluant along the sheet.

It is therefore necessary to insert the sheet of chromatography paper so that the fibers run in the same direction as the migration of the solvent.

It is easy to determine the orientation of the fibers in chromatography paper by tearing a piece in one direction and then in the other direction. The paper will tear preferentially in the direction corresponding to the orientation of the fibers.

4 Maintenance

Cleaning

The parts should be cleaned in soapy water. Never use organic solvents, as they may damage the parts.

To make cleaning the vessel easier and to avoid the formation of mineral deposits after drying, use demineralised water for the final rinse.

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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Service gratuit * :
+ 33 (0)2 32 29 40 50

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Rubrique FAQ

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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