# **Ex.A.O**

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## ***Data logging***

*REDY Wireless Sensor*

Réf :  
**488 032**

Français – p 1

English – p 14

Version : 3207

***Capteur REDY Poulie***  
***REDY Stepped Pulley Sensor***

## Table des matières

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Table of Contents

1. Content .....	15
2. Description.....	15
2.1 REDY sensor .....	15
2.2 Pulley .....	16
3. First Use .....	16
4. Starting Up.....	16
5. Using the Sensor in Connected Mode to a Windows Computer.....	16
5.1 Installing the software Redy for Windows .....	17
5.2 Connecting the Sensor to the Software.....	17
5.3 Mechanical calibration of the pulley.....	18
5.4 Carrying Out an Acquisition Depending on Time .....	18
6. Using the Sensor in Bluetooth Mode on Tablet .....	22
6.1 Installing the application Redy for Tablet.....	22
6.2 Connecting the Sensor to the Tablet via Bluetooth .....	22
6.3 Setting the Acquisition Depending on Time .....	24
6.3.1 Setting the Intensity Measurement.....	24
6.3.2 Setting the Time Factor .....	24
7. Technical Specifications.....	26

## REDY Stepped Pulley Sensor


Ref : 488 032

### 1. Content


- 1 angular speed / angle / distance Sensor REDY
- 1 Stepped Electronic Pulley attached to the Redy Sensor
- 1 type C USB cable
- 1 protective Cover

### 2. Description

#### 2.1 REDY sensor

## 2.2 Pulley

## 3. First Use

When receiving your sensor, and before the first use, it is recommended to fully charge the battery.

Use the USB-C cable supplied to connect the sensor to a USB port of a computer, or a power outlet via a universal charger (like the USB mains charger ref. 805819, not supplied)

Before each practical work, plan to fully charge the device.  
The battery life varies depending on the conditions of use and the number of measurements and is designed to ensure a full practical session.

## 4. Starting Up

Press the Power button.  
The device turns on. The value measured appears on the screen.  
Your device is ready for use.

## 5. Using the Sensor in Connected Mode to a Windows Computer

### 5.1 Installing the software Redy for Windows


The software Redy for windows is available for downloading on the Jeulin website.

### 5.2 Connecting the Sensor to the Software

- Connect the sensor to a USB port of your computer.
- Turn on the sensor by pressing on the position I, of the Power button.
- Launch the software Redy for Windows by double-clicking on the icon present on the desktop of the computer.

The sensor is automatically detected. Its symbol appears in the software list of sensors :

## REDY Stepped Pulley Sensor

Ref : 488 032

### 5.3 Mechanical calibration of the pulley

Before starting the measurements, it is important to perform the mechanical zero adjustment.


Place the pulley on its support.

Attach all accessories, wait until there are no more movements.

Slightly loosen the fixing and orient the pulley so that the index finger is as far as possible from the zero of the serigraphy.


Retighten the fastener.

Once the index is opposite zero and the pulley is no longer moving, press the calibration button on the back of the pulley.


### 5.4 Carrying Out an Acquisition Depending on Time


To carry out an acquisition with the sensor, drag and drop the icon assigned to the desired measurement (angle, angular speed or distance), in ordinate using the mouse :

All the sensor relative indications then appear under the graph:


By clicking on the tab Measure, the value measured by the sensor appears on the screen:

# REDY Stepped Pulley Sensor


Ref : 488 032

On the Redy software, after entering the ordinates, click on Setting, then on Set.  
 Note that calibration must be done before making angle or distance measurements. Do not forget to select the measurement to be calibrated in the y-axis.


Calibration is complete.

It remains to set the time, by dragging and dropping the timer icon in abscissa in the graph, as shown below :

All the indications relative to setting the time acquisition then appear under the graph :

Hence, the user can, if he likes, modify the acquisition duration and/or the number of points.

Once the configurations are done, start the acquisition by pressing on the button « Start » (see above).

## 6. Using the Sensor in Bluetooth Mode on Tablet

### 6.1 Installing the application Redy for Tablet


The application Redy is free and available on Google Play or the App Store.

### 6.2 Connecting the Sensor to the Tablet via Bluetooth

- Turn on the sensor by pressing on the position I, on the Power button.
- On your tablet, launch the application Redy, by pressing on the icon :

- The application is launched, you are asked to connect a sensor :

# REDY Stepped Pulley Sensor

Ref : 488 032

- Press on yes.

The sensors management space opens automatically in the application.

The sensor appears in the sensors list.

In order to make sure that the desired sensor is detected, you can verify that the address number written on the sensor screen, and that appearing in the application sensors list, match.

This address number is unique for each sensor.

*In the example above, it is the "000 - Rotary motion".*

- Start the connection by pressing on the desired sensor in the list. A temporary icon appears, then the application confirms that the sensor is connected when the indications below show on the screen :

Note the apparition of the Bluetooth logo next  to the sensor name, as well as the message « Connected » in the end.

- Go back to the main screen by pressing on the back arrow :

## 6.3 Setting the Acquisition Depending on Time

### 6.3.1 Setting the Intensity Measurement


The sensor activation, detailed in the previous chapter, selects the sensor automatically for an acquisition. Therefore, no further additional configuration action needs to be taken.

### 6.3.2 Setting the Time Factor


Click on the icon Setting of the main screen :

Then in the window Acquisition, select Time :


The user can, if he likes, modify the duration and the time unit.

Go back to the main screen by pressing on the back arrow :

Then, start the acquisition by pressing on the Start button :

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

**Free service\***

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux