**Ex.A.O**

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***Data logging***

*REDY Wireless Sensor*

Réf :  
**488 012**

Français – p 1

English – P 12

Version : 4204

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**REDY Ohmmeter Sensor**

## Table des matières

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
## Table des matières

1. Content .....	13
2. Description.....	13
3. First Use .....	14
4. Starting Up .....	14
5. Using the Sensor in Connected Mode to a Windows Computer .....	14
5.1 Installing the software Redy for Windows .....	14
5.2 Connecting the Sensor to the Software.....	14
5.3 Carrying Out an Acquisition versus Time.....	15
6. Using the Sensor in Bluetooth Mode on Tablet .....	17
6.1 Installing the application Redy for Tablet.....	17
6.2 Connecting the Sensor to the Tablet via Bluetooth .....	17
6.3 Display of the measurement unit .....	18
6.4 Setting the Acquisition versus Time .....	20
6.4.1 Setting the measurement.....	20
6.4.2 Setting the Time parameters.....	20
7. Technical Specifications.....	22

## 1. Content

- 1 Ohmmeter sensor REDY
- 1 type C USB cable
- 1 protective cover

## 2. Description

### 3. First Use

When receiving your sensor, and before the first use, it is recommended to fully charge the battery.

Use the USB-C cable supplied to connect the sensor to a USB port of a computer, or a power outlet via a universal charger (like the USB mains charger ref. 805819, not supplied)

Before each practical work, plan to fully charge the device.

The battery life varies depending on the conditions of use and the number of measurements and is designed to ensure a full practical session.

### 4. Starting Up

Press the Power button.

The device turns on. The value measured appears on the screen.

Your device is ready to use.

## 5. Using the Sensor in Connected Mode to a Windows Computer

### 5.1 Installing the software Redy for Windows


The software Redy for windows is available for downloading on the Jeulin website.

### 5.2 Connecting the Sensor to the Software

- Connect the sensor to a USB port of your computer.
- Turn on the sensor by pressing on the position I, of the Power button.
- Launch the software Redy for Windows by double-clicking on the icon present on the desktop of the computer.


The sensor is automatically detected. Its symbol appears in the software list of sensors:


### 5.3 Carrying Out an Acquisition Depending on Time

To carry out an acquisition with the sensor, drag and drop the sensor icon, in ordinate using the mouse:


All the indications relative to the sensor then appear under the graph:


It remains to set the time, by dragging and dropping the timer icon in abscissa in the graph, as shown below:

All the indications relative to setting the time acquisition then appear under the graph:

Once the configurations are done, start the acquisition by pressing on the button « Start » (see above).


## 6. Using the Sensor in Bluetooth Mode on Tablet

### 6.1 Installing the application Redy for Tablet


The application Redy is free and available on Google Play or the App Store.

### 6.2 Connecting the Sensor to the Tablet via Bluetooth

- Turn on the sensor by pressing on the position I, on the Power button.
- On your tablet, launch the application Redy, by tapping on the icon:


- The application is launched, you are asked to connect a sensor:

- Press on yes.

The sensors management space opens automatically in the application.


The sensor appears in the sensors list.

In order to make sure that the desired sensor is detected, you can verify that the address number written on the sensor screen, and that appearing in the application sensors list, match.

This address number is unique for each sensor.

*In the example above, it is the 000 – Ohmmeter.*

- Start the connection by tapping on the desired sensor in the list. A temporary icon appears, then the application confirms that the sensor is connected when the indications below show on the screen:

Note the apparition of the Bluetooth logo  next to the sensor name, as well as the message « Connected » at the end.

### 6.3 Display of the measurement unit

Click on the sensor settings button in the window below:

Choose mode and language, if need be :

 Display

$R (\Omega)$   
512


Mode

10k $\Omega$  : R (k $\Omega$ )

Language

Select

 Display

$R (\Omega)$   
2

Mode

10k $\Omega$  : R (k $\Omega$ )  
100k $\Omega$  : R (k $\Omega$ )  
ANNULER

Mode

10k $\Omega$  : R (k $\Omega$ )

Language

Select

Return to the main screen by tapping the back arrow twice:

 Display

$u (mV)$   
-0.20


Mode

Voltage : u (+/-1 mV)

Language

English

The value read by the sensor appears in the multimeter screen, placed by default in the bottom right portion of the screen:


## 6.4 Setting the Acquisition versus Time

### 6.4.1 Setting the measurement


The sensor activation, detailed in the previous chapter, selects the sensor automatically for an acquisition. Therefore, no further additional configuration action needs to be taken.

### 6.4.2 Setting the Time parameters

Click on the icon Setting of the main screen:

Then in the window Acquisition, select Time:


The user has to modify the duration and the time unit.

A screenshot of the 'Acquisition : Time' configuration screen. The screen has a dark red header with a back arrow and the title 'Acquisition : Time'. Below the header is a blue bar with 'Duration' and 'Synchronization' labels. The 'Duration' field is set to '10' and the 'Unit' is set to 'Second'. There is a 'Continuous acquisition' toggle switch at the bottom, which is currently turned off. A small text note below the toggle says: 'Allows the acquisition to loop to infinity, this option is useful for short measurements (measurement of a sound for example).'.

Go back to the main screen by tapping on the back arrow:

A screenshot of the 'Acquisition : Time' configuration screen, identical to the one above. A red circle highlights the back arrow in the top left corner of the dark red header bar.

Once the configurations are done, start the acquisition by pressing on the button « Start ».

## 7. Technical Specifications

- Backlit graphic display.
- Li-ion rechargeable battery.
- USB-C connection.
- Bluetooth 4.0 connection

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediately to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

### Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux