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Spectrophotométrie

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
**0 825 563 563\***

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## 1. Presentation

This USB/Bluetooth spectrophotometer has been developed for measurements :

In physics :

- Produce and analyze emission spectra from internal or external light sources (using the fiber optic accessory).
- Discover the notions of absorbance/transmittance and complementary colors, thanks to a slot for colored filters.

In chemistry :

- Measuring absorbance and transmittance of solutions
- Measurements as a function of concentration or time
- Measurement and analysis of emission spectra of luminescent solutions not subjected to an excitation wavelength (chemiluminescence), thanks to the possibility of extinguishing the integrated light source.

In SVT :

- Measurement and analysis of emission spectra of fluorescent compounds subjected to an excitation wavelength (chlorophyll).

### 1.1 Composition

#### **Spectrophotometer only Réf.705414**

Spectrophotometer  
 Tank/filter cover  
 USB Cable  
 Power supply

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Spectrophotometer  
 Tank/filter cover  
 USB Cable  
 Power supply  
 Optical fiber

## 1.2 Technical data

<b>Measurements modes :</b>	Transmittance Absorbance Fluorescence/Chemiluminescence
<b>Resolution :</b>	2 - 3 nm
<b>Samples :</b>	10x10mm Tank ou tube 50x50mm Slide
<b>Spectral range :</b>	360 – 1000 nm
<b>Light sources :</b>	White LED + selectable colors (via software) Lateral excitation LEDs 395-415 nm et 500-540 nm
<b>Connection :</b>	USB Bluetooth
<b>Software :</b>	Intégrated (USB mode) – Plug & Play
<b>Application :</b>	Download (Bluetooth mode)
<b>Power supply :</b>	Built-in battery + charger unit included
<b>Python compatibility :</b>	Yes
<b>Software compatibility :</b>	Windows 7 and above, Mac, Android, iOS
<b>Accessories included :</b>	Tank/filter cover USB cable Power supply 12V-1A
<b>Battery :</b>	Use : 8 hours Maximum charge time : 4 hours 190 x 120 x 60 mm

## 2. Start


For maximum ease of use, this device features "plug and play" software that launches automatically as soon as your device is connected to the USB port.

It is also possible to install the software on a PC (see 2.2).

This will enable you to use the software without having to plug in the spectrophotometer (for example, to process curves during an exercise, or to write a report...).

Note: similar navigation and software functions via USB and Bluetooth.


Help is also integrated into the software on the PC, and can be accessed via the Help tab:


## 2.1 Lunch in « Plug & Play » mode

- Switch on your computer.
- Connect the device to the computer using USB cable included.
- The following message will appear, click on **Run spectro.exe** to lunch the software.


## 2.2 Lunch with software installation

- Connect the spectrophotometer to the computer, and it twill appear in **Computer** or **This PC** , a **Spectro** drive:

- Right-click on the drive, choose **Open**, then in the installation folder, double-click on **Spectro\_EN.msi**


- The software installation starts, juste follow the instructions to complete the installation.

## 2.3 Internal memory

Files (curves, protocols, photos, etc.) can be saved on your computer, on the school network or in the spectrophotometer itself.  
 In fact, the spectrophotometer has an internal memory that enables you to save your documents in the instrument.


This internal memory can be accessed in **Computer** or **This PC**, under the name **My\_Spectro**:


When you open this drive, you can drag-and-drop (or copy-and-paste) various documents, just as you would on a USB stick.

It is also possible to save curves at the end of acquisition.


From the spectrophotometer software, simply select **File** then **Save as** :

Then select the location indicated above (**This PC > My\_Spectro**)


To open a file saved in the spectrophotometer, simply select **File** then **Open**.

## 3. Use with integrated software (USB)

### 3.1 Emission tab

Emission mode measures the emission spectrum of light sources.

- **Integration time**

It allows you to modify the gain during acquisition. Increasing the integration time increases the intensity, but moves the zero point.

Caution: increasing the integration time also increases the noise intensity.

The checkbox next to the value lets you manage integration time automatically or manually. In manual mode (box unchecked), you can enter the desired value directly in the display field.

- **Black measurement**

This function allows you to measure a reference curve of the device.

Click on the "Black measurement" button, and follow the instructions. The operation may take a few seconds.

- **Freeze curve [Pause]**

The spectrum is acquired in real time. It is necessary to suspend acquisition (and therefore freeze the curve) to be able to use all the software's functions to evaluate the results or save your curve.

To do this, use the Play/Pause button:

To freeze the curve

To resume the acquisition

**Please note:** when the curve is frozen, it is not saved, so you need to keep it (see below).

- **Keeping a curve**

Click on

Name the curve and it will be saved.

#### Notes :

- The curve acquired in real time is called **I( $\lambda$ )**, and can be displayed or hidden like the other curves.
- Don't forget to click **Play** again after saving the curve.
- If you click directly on **Keep**, the curve is paused for the duration of the curve recording.

#### • Curve smoothing

The **Average** and **Smooth** options activate mathematical processing functions to reduce irregularities in the curve:


- **Average** to eliminate noise, adjustable as *No - Low - Medium - High*
- **Smooth**, has to be limited so as not to lose too much resolution, adjustable in *No - Low - Medium - High*.

Please note: The smoothing parameters selected in the emission tab are set by default for the other tabs.

#### • Obtaining a color spectrum

The spectrum can be displayed as an envelope or as a colored spectrum. To do this, use the following function.

☒ Color the area under the curve

## 3.2 Absorbance tab

The absorbance mode measures the light intensity received after passing through your sample (cell or filter).

Before measuring absorbance or transmittance, you need to perform a calibration. To do this, insert a cuvette filled with solvent and click on **Calibrate**, the instrument acquires the reference spectrum. After calibration, the black level is measured. This determines the level of noise (optical, electronic...) to be subtracted from the curve during measurement.

Insert the colored filter or the cuvette containing your colored solution. The corresponding spectrum is displayed in real time in the graphics area.

Data  
 Emission  
**Absorbance / Transmittance**  
 Fluorescence

### Absorbance / Transmittance

Pause
 Save

Integration time 19.4 ms ☒ Automatic

Calibrate

**Display**


☐ Intensity
 ☒ Absorbance
 ☐ Transmittance

**Wavelength**

nm

Concentration
 kinetic

Example of a curve obtained for a potassium permanganate solution :

**Note:** the procedure for freezing, storing and displaying curves is the same as for the Emission section.

The software can then be used to check Beer's law directly, without the need for an external spreadsheet.

This part follows on from the previous paragraph, in which calibration was carried out.

① Select the working wavelength using the slider, or by entering it manually in the box provided.

Now choose the type of experiment to be performed:

② Concentration: for Beer-Lambert-type measurements,

③ Kinetics: for time-dependent measurements

### 3.2.1 Absorbance tab > Concentration

Indicate the quantity measured on the x-axis (Volume, concentration...etc) and its unit, generally the concentration of the solutions analyzed.


For each sample, indicate the value measured on the x-axis, then click on **Validation**.

After each sample has been acquired, a point is added to the graph showing the absorbance of the solution as a function of its concentration.

This curve directly illustrates Beer-Lambert's law  $A = f(C)$ . The calibration line can be drawn automatically by clicking on **Draw calibration line**.

This line is calculated by linear regression on all acquired points.

Place the solution of unknown concentration in the cuvette holder. The software measures the absorbance and deduces the concentration. Its values are displayed in the **Unknown solution** box.

### 3.2.2 Absorbance tab > Kinetic

For the study of slow reactions causing a change in colored solution, such as the oxidation of iodide ions, the spectrophotometer can be used to determine the kinetics of this reaction.

To do this, select **Kinetics**. the window shown opposite appears.

The kinetics can be studied at two different wavelengths. If the study is at a single wavelength - a single value is expected.

Select the time interval between two acquisitions, then start the acquisition.

You can then stop acquisition by clicking on **Stop**.

In the case of automatic stop, you can pre-select the acquisition duration. It will automatically stop after the set time (in seconds).

Stop

Time between two measurements

1 s

### 3.3 Fluorescence tab

The fluorescence mode measures the fluorescence emitted by your sample following its excitation by a selected wavelength.

Two wavelength ranges are available. It is also possible to emit no excitation and still collect the spontaneous fluorescence emission from the sample.

Before any fluorescence measurement, you need to perform a calibration. To do this, insert a cuvette filled with solvent and click on **Calibrate**, the instrument acquires the reference spectrum. Blackness measurement follows this calibration. This step determines the level of noise (optical, electronic...) to be subtracted from the curve during measurement.

### 3.3.1 Fluorescence tab > Concentration

Indicate the quantity measured on the x-axis (Volume, concentration...etc) and its unit, generally the concentration of the solutions analyzed.


For each sample, indicate the value measured on the x-axis, then click on **Validation**.

After each sample has been acquired, a point is added to the graph showing the fluorescence of the solution as a function of its concentration.

The calibration line can be drawn automatically by clicking on **Draw calibration line**.

This line is calculated by linear regression on all acquired points.

Place the solution of unknown concentration in the cuvette holder. The software measures the absorbance and deduces the concentration. Its values are displayed in the **Unknown solution box**.

### 3.3.2 Fluorescence tab > Kinetic


For the study of slow reactions causing a change in colored solution, the spectrophotometer can be used to determine the kinetics of this reaction

To do this, select **Kinetics**. the window shown opposite appears

Select the time interval between two acquisitions, then start the acquisition.



You can then stop acquisition by clicking on **Stop**.

In the case of automatic stop, you can pre-select the acquisition duration. It will automatically stop after the set time (in seconds).

## 4. Use with downloadable application (Bluetooth)

You can download the application from the Google PlayStore or App Store, or scan the QRCode on the spectrophotometer :

The applications and software are designed so that the menus match each other. As a result, it is possible to rely on the previous elements of this manual concerning the software to carry out measurements with ease.

## 5. After-Sales Service

This material is under a two years warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**+33 (0)2 32 29 40 50**