Prélèvement

Aspiration

Ref :
703 323
703 326
703 327
703 381

Français – p 1

English – p 5

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1 Introduction

These micropipettes are for general use, for aspiration and dispensing precise volumes of liquids. They are continuously adjustable. They work on the air displacement principle (the volume of air between the plunger and the liquid). The volume of liquid to be sampled can be adjusted with an adjustment dial, and it is numerically displayed and it is easy to read, thereby reducing the risk of errors. The proper disposable, polypropylene cone tips must be used with these devices. The size and the volume of these cone tips depend on the micropipette and the range of volumes to be aspirated.





Model	Item no.	Volumes
H20	703 323	2 - 20 μ L
H100	703 326	10 - 100 μ L
H200	703 381	20 - 200 μ L
H1000	703 327	100 - 1000 μ L

2 Technical data

- Adjustable volume
- Numerical display of the selected volume
- Cone tip ejection system
- Separate cone tip ejection button located on the side

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4 Operation

4.1 Adjusting the volume

The volume desired is set by turning the pushbutton clockwise to lower the volume and counter clockwise to increase it.

For each of the pipettes:

- never attempt to set volumes outside the range specified in the above table. Doing so will affect the precision of the mechanism;
- never attempt to aspirate volumes outside the range specified in the above table, as the volumes thus aspirated will be inaccurate. If you aspirate too large a volume, you risk causing the liquid to penetrate the device and permanently contaminating your material.

4.2 Ejecting the cone tips

To reduce the risks of contamination, each pipette is equipped with a cone tip ejection system. The system consists of a pushbutton located on the side. Point the pipette towards the waste receptacle. Press the button on the side with your thumb; the cone tip will be ejected safely.

4.3 Pipetting

Warning:

- For greater precision, push and release the pushbutton gently, especially with highly viscous liquids. Never release the pushbutton suddenly.
- Optimum use temperature: 20°C ± 5°C
- Firmly attach a clean, suitable cone tip to the collar of the pipette.
- Set the desired volume to be pipetted.
- Gently press the pushbutton down to the **first** stop.
- Immerse the cone tip of the pipette vertically in the solution to be pipetted (to a depth of about 2 to 4 mm).
- **Slowly** let the pushbutton return to its original position; the liquid will be aspirated slowly.
- After 1 to 2 seconds, remove the pipette from the liquid.

Aspiration
Adjustable volume micropipettes
Ref :
703 323
703 326
703 327
703

4.4 Dispensing the liquid

- Place the cone tip of the pipette in the container to receive the liquid being dispensed.
- Place the cone tip near the wall of the container.
- Press the pushbutton down to the first stop, and after 1 second press it to the second notch to empty all of the liquid contained in the cone tip.
- Remove the pipette from the liquid after 1-2 seconds.
- **Slowly** let the pushbutton return to its original position.

5 Maintenance

The device cannot be disassembled for maintenance.
Non-autoclavable.
If possible, store vertically. Store in a dry place.

6 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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