
Overhead Stirrer SS10 & SS20

Instruction Manual Version 1.4


Figure 1

Thank you for purchasing this Stuart product. To get the best performance from the equipment, and for your personal safety, please read these instructions carefully before use.

This equipment is designed to operate under the following conditions:-

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuation not exceeding 10%
- ❖ Overvoltage category II IEC60364-4-443
- ❖ Pollution degree 2 IEC664

If the equipment is not used in the manner described in this manual the protection provided by the equipment may be impaired.

WARNING! This Overhead Stirrer is not suitable for use in inflammable or explosive atmospheres. It should not be used to stir any inflammable liquids unless they are in a suitable sealed vessel where there is no possibility of gases or vapours coming into contact with the Overhead Stirrer body.

Electrical Installation

THIS EQUIPMENT MUST BE EARTHED

Before connection please ensure that the mains power supply corresponds to that shown on the rating plate.

All models are supplied with a mains lead fitted with an IEC plug to fit the IEC socket on the short flying lead emerging from the Stirrer mounting bar. This lead should be connected to the flying lead BEFORE connection to the mains supply. Should the mains lead require replacement a cable of 1mm² of harmonised code H05W-F connected to an IEC 320 plug should be used.

Most Stuart equipment is supplied with the mains lead fitted with a moulded plug for connection to the electricity supply. If this is not the case it should be fitted with a plug approved for the country in which it is to be used, noting that the mains lead has cables coloured according to the following code:

Blue-	Neutral
Brown-	Live
Green/yellow-	Earth

If required the plug should be fitted with a 10A fuse to protect the cable. All models are protected by two fuses located in the back panel. These conform to IEC127 and are rated as follows:

Model	Power	Fuse
SS10	50W	T1A
SS20	80W	T1A

IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN

General Description

The Stuart Overhead Stirrers are designed for effective mixing of large volumes or viscous liquids. All models have an easy to use keyless chuck, an enclosed mains cable designed to be kept away from moving parts, and overload protection for maximum safety.

SS10 General Purpose Stirrer: is designed for use in aqueous and low viscosity liquids to a maximum volume of 15 litres.

SS20 High Performance Stirrer: is designed for use in low to medium viscosity liquids to a maximum volume of 25 litres. It has a hollow spindle to enable long stirrer shafts to be used or for raising the paddle out of vessels easily.

Viscosity	mPa-s	Example
Very low	0 - 100	Water to thin oil
Low	100 - 1000	Thin to thick oil
Medium	1000 - 10,000	Thick oil
High	10,000 - 20,000	Thick oils to honey

Safety Advice Before Use

- ❖ Never operate the Stirrer without first clamping the mounting rod securely. It is NOT intended to be hand held. In use the stability of the Overhead Stirrer is dependant on the security of the mounting assembly. The operator should ensure that the mounting is adequate for the intended application. (See the Accessories Section for stand and stirrer paddles available).
- ❖ Clamp the stirrer shaft securely in the chuck and lock before operating the Stirrer.
- ❖ Do not operate model SS20 or with a stirrer shaft protruding through the top of the Overhead Stirrer without first removing the top grommet.
- ❖ Do not mix flammable liquids or use the equipment in hazardous atmospheres.
- ❖ In case of mains interrupting, a fault or mechanical failure, the unit will continue to operate on removal of fault.
- ❖ Mechanical energy can lead to breakage of glass vessels. Use with care.

WARNING! The Overhead Stirrer will automatically restart after a power failure. It may also restart in the event of a fault or mechanical interruption (stall).

Operation

Fixing a Stirrer Paddle

The chuck is designed to grip a stirrer shaft up to a maximum of 13mm diameter. Loosen the jaws of the chuck by turning clockwise. Place a stirrer shaft in the chuck and tighten by turning anti clockwise. The chuck only requires tightening by hand.

The chuck protective cover is designed for use with stirrer shafts up to 8mm diameter only. It grips the top of the shaft preventing ingress of corrosive materials. If a stirrer shaft with a diameter greater than 8mm is used, remove the protective cover before fixing the shaft in the chuck.

Hollow Spindle (SS20 only)

Model SS20 has a hollow spindle enabling a stirrer shaft of diameter up to 8mm to be passed through the Stirrer body and adjusted from the top. Open the chuck and remove the top grommet. Feed the shaft through the chuck and into the Stirrer until it protrudes through the top. Tighten the chuck as described above.

Using the Stirrer (SS10 & SS20)

Clamp the Overhead Stirrer securely to a stand and fix a stirrer paddle in the chuck (as described above). Check that the Stirrer is switched OFF before inserting the mains cable into the short flying lead and connecting to the mains electricity supply. When the supply is available from the mains, a green light will illuminate on the front of the Stirrer.

Selecting the Stirring Speed (SS10 & SS20)

Switch the Stirrer ON by the mains on / off switch. A yellow light will illuminate on the front of the Stirrer. Rotate the speed control knob to select the required speed of rotation. A soft start facility is built into the Stirrer so there will be a slight delay between selection and rotation.

Overload Protection (SS10 & SS20)

If the load on the Overhead Stirrer is too great a red warning light on the front of the Stirrer will illuminate and the power to the motor will be limited to protect the Stirrer from damage. The Stirrer will continue to run at the limited speed until the load is reduced.

If there is a danger of overheating the power to the motor will be further reduced. If the danger of overheating becomes severe the Stirrer will stop. The unit must be switched OFF by the mains on / off switch and allowed to cool before it can be used again.

Maintenance & Servicing

WARNING: Ensure the unit is disconnected from the mains electricity supply before attempting maintenance or servicing.

Periodically clean the instrument using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents.

Any repairs or replacement of parts MUST be undertaken by suitably qualified personnel.

Spares

The following spares and accessories are available from your laboratory dealer:

Description	Catalogue Number
Fuse, T1A	ESS20C1(S)
Top case grommet	USS20C5
Chuck protective cover	USS2056
For complete paddle, order rod plus head(s)	
Paddle rod, 350mm	SS10/5
Paddle rod, 550mm	SS10/6
Small paddle head	SS10/10
Large paddle head	SS10/11
Large propeller head (2 blade)	SS10/12
Small propeller head (4 blade)	SS10/13
Anchor paddle head	SS10/14
Stand, General Purpose (inc. base, rod & bosshead)	SS10/1
Stand, Heavy Duty (inc. base, rod & bosshead)	SS10/2

For a comprehensive list of parts and a Service Manual required by service engineers conducting internal repairs, please contact the Technical Service Department of Cole-Parmer quoting both the model and serial number.

Only spare parts supplied or specified by Cole-Parmer or its agents should be used. Fitting of non-approved parts may affect the

performance and safety features designed into the instrument.

If in any doubt, please contact the Technical Dept. of Cole-Parmer or the point of sale.

Cole-Parmer Ltd.
Beacon Road,
Stone, Staffordshire,
ST15 0SA, United Kingdom
Tel: +44 (0)1785 812121
Email: cpservice@coleparmer.com
Web: www.stuart-equipment.com

Warranty

Cole-Parmer Ltd. warrants this instrument to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of three (3) years. In the event of a justified claim, Cole-Parmer will replace any defective component or replace the unit free of charge.

This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by installation, adaptation, modification, fitting of non-approved parts or repair by unauthorized personnel.

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Cole-Parmer Ltd.
Beacon Road,
Stone, Staffordshire,
ST15 0SA, United Kingdom
Tel: +44 (0)1785 812121
Email: cpSERVICE@coleparmer.com
Web: www.stuart-equipment.com

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Notes

Notes

This product meets the applicable EC harmonized standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in its vicinity will meet these standards

and so we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, contact the manufacturer.

EU Declaration of Conformity

Product	Laboratory Equipment	File Number	P225
Manufacturer	Cole-Parmer Ltd Beacon Road Stone, Staffordshire ST15 0SA United Kingdom		

This declaration of conformity is issued under the sole responsibility of the manufacturer

Object of Declaration Overhead Stirrers
(reference the attached list of catalogue numbers)

The object of the declaration described above is in conformity with the relevant Union Harmonisation Legislation:

Low Voltage Directive	2006/95/EC	(until 19 April 2016)
	2014/35/EU	(from 20 April 2016)
EMC Directive	2004/108/EC	
RoHS Directive	2011/65/EC	

References to the relevant harmonised standards used or references to the other technical specifications in relation to which conformity is declared:

IEC/EN 61010-1:2010	Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1: General requirements.
IEC/EN 61010-2-051:2015	Particular requirements for laboratory equipment for mixing and stirring.
IEC/EN 61326-1:2006	Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 1: General requirements (Class A).

Signed for and on behalf of the above manufacturer

Additional Information	Year of CE Marking: 2001
Place of Issue	Stone, Staffordshire, UK
Date of Issue	June 2009. Revised 06 February 2017
Authorised Representative	Carl Warren
Title	Technical Manager
Signature	

INSPECTION REPORT

MODEL _____

ELECTRICAL SAFETY

1. Earth continuity

2. Insulation

3. Flash test

FUNCTIONAL

1. Indicators


2. Speed action

3. Visual acceptance

QUALITY CONTROL INSPECTOR

Cole-Parmer Ltd.

Beacon Road,
Stone,

Staffordshire,
ST15 0SA,

United Kingdom

Tel: +44 (0)1785 812121

Email: cpinfo@coleparmer.com

Web: www.stuart-equipment.com