# Mécanique

Ultrasons

# *Mechanics*

*Ultrasounds*

Réf :  
222 094

Français – p 1

English – p 5

Version : 2206

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## 1. Description

This note describes the functionalities of this transmitter telemeter Moduson®3 and the compatibilities with the Moduson®3 range.

The transmitter telemeter (ref 222 077) powered by 12V power supply part. 281 612.

The frequency of the transmitter can be adjusted to about 40 kHz.

The transmission of ultrasonic waves is triggered by an On/Off switch.

A reversing switch helps obtain either a continuous transmission or a transmission in bursts. A second reversing switch helps select the duration of the bursts:

- Long bursts: duration of 6 ms, with a gap of 4 ms between bursts
- Short bursts: duration of 1 ms, with a gap of 4 ms between bursts

This function is used to measure distances from 0 to 30 cm using the sonar and scan techniques. Actually, this technique requires that the echo reflected by the object does not get mixed up with the burst emitted, which requires the use of short bursts.

## 2. Ideal configuration


The ideal configuration for the measurements is the following:

- 1 transmitter telemeter Moduson®3 part. 222 094
- 1 external power supply 12V 600mA par. 281 612
- 2 receivers Moduson®2 part. 222 078 (for measurements in the air only) or Moduson®3 part. 222 096 (for measurements in the air or in a fluid)
- 1 tank with a support, in case of liquid measurements, contained in the pack part. 222 097
- 1 ExAO console or an oscilloscope


## 2.1 Use with receivers Moduson®3 and the Foxy console

In this part, we will explain how to use the transmitter telemeter Moduson®3 with one or two receivers Moduson®.


The transducers of the new versions of Moduson® allow a use in the fluids and in the air.

That is why if there is a change in the environment, it is necessary to adjust the emission gain (therefore the reception) of the signal in order to obtain signals which are not deformed in case of a significant gain.

Example of a set-up – Measurement of the speed sound in water

## 2.2 Use with receivers Moduson®2 and the Foxy console

In the case of using receivers Moduson®2, no connection besides the power supply and the BNC/Banana plug has to be done.

In fact, it will be necessary to connect the receivers Moduson®2 directly on the Foxy console, because they have connections such as safety banana, of diameter 4mm.

### 3. Instructions of use

1. In order to guarantee optimum signals, it is necessary to ensure that the transmitter Moduson®3's mass and the Foxy console's mass are connected.
2. The transmitter telemeter has a self-resetting fuse, in case the connection creates a short circuit, the device puts itself in safety and does not produce any signal. In this case, turn the device's power off, eliminate the short-circuit, wait a few seconds and switch it back on.
3. Adjusting the gain to overcome a distorted signal in the air.


### 4. Examples of possible experiments

- 1 Highlight the vibrational nature of an ultrasound, measure its period and frequency in the air or in a fluid
- 2 Measure an ultrasound in the air or in a fluid
- 3 Observe the existence of points that vibrate in phase
- 4 Verify the relation  $\lambda = v/N$
- 5 Measure a distance by ultrasounds in the air or in a fluid

### 5. Examples of results

#### 5.1 Ultrasounds speed in the air

Results:

Distance between the 2 receivers: 14 cm  
 Interval between the 2 receivers: 0.4113 ms  
 Reminder of the theoretical speed of sound in air: 340 m.s<sup>-1</sup>

$$\Delta d = 14 \text{ cm}$$

$$\Delta t = 0.0004113 \text{ s}$$


$$v = \frac{\Delta d}{\Delta t} = \frac{0.14}{0.0004113}$$

$$v = 340 \text{ m.s}^{-1}$$


## 5.2 Speed of the ultrasounds in water

Positioning of the equipment:

Results:

Distance between the 2 receivers: 18 cm

Interval between the 2 receivers: 121.9  $\mu$ s

Reminder of the theoretical speed of sound in water: 1460 m.s<sup>-1</sup>

$$\begin{aligned}\Delta d &= 18 \text{ cm} \\ \Delta t &= 0.0001219 \text{ s} \\ v &= \frac{\Delta d}{\Delta t} = \frac{0.1800}{0.0001219}\end{aligned}$$

$$v = 1476 \text{ m.s}^{-1}$$

## 6. After-sales service

The device is under a 2-year guarantee.

For any repairs, contact the **Technical Support** on **+33 (0)2 32 29 40 23**.

The equipment should be returned in our workshops and for any repairs or spares parts, please contact:

**JEULIN – TECHNICAL SUPPORT**  
 468 rue Jacques Monod  
 CS 21900  
 27019 EVREUX CEDEX FRANCE  
**+33 (0)2 32 29 40 23**