Spectrophotométrie

Kit didactique

Spectrophotometry

Didactic kit

Réf :
105 479

Français – p 1

English – p 9

Version : 4206

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+33 (0)9 69 32 02 10

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1. Introduction

This kit has been designed as an introduction to the observation and measurement of fluorescence in spectrophotometry.

The compounds in the kit have been selected for their specific fluorescence characteristics and behaviors.

This kit has been developed for use with the Jeulin USB/Bluetooth spectrophotometer, but can be adapted to other models.


The dosages indicated make it possible to carry out the initiation manipulations, but there is enough material to do other types of TP (concentration range etc.)

Some compounds are very dirty, wearing classic PPE (goggles, gloves, blouse) is recommended.

2. Composition

This kit includes :

- 1 x microtube containing approx. 0.2 g Riboflavin
- 1 x microtube containing approx. 0.1 g Quinine
- 1 x microtube containing approx. 0.8 g Fluorescein
- 1 x microtube containing approx. 0.8 g Eosin
- 1 x 60 mL containing approx. 0.2 g Rhodamine solution in 60mL distilled water

3. Experiment

3.1 Riboflavin

3.1.1 Introduction

Riboflavin, or vitamin B2, is essential for energy metabolism and the health of the skin, eyes and nervous system. Found in energy drinks, it helps convert nutrients into energy. It is an essential component of many coenzymes required for cellular respiration.


3.1.2 Use

The microtube contains approx. 0.2 g Riboflavin, dose for one liter.
To observe fluorescence and set up a measuring range, accurately weigh and dilute the entire microtube in 1L of distilled water.

This solution will be used to study the compound's fluorescence response.
An energy drink dosing session can also be carried out, by making a standard range from the stock solution.

Solution storage: Store for several days in a dark place.
Long-term storage will result in a drop in fluorescence intensity.

3.1.3 Results

Excitation

Absorbance maximum 448 nm
Spectrophotometer setting 405 nm

Emission

Fluorescence maximum 530 nm

3.1.4 Data

Name : Riboflavin 98%
CAS : 83-88-5
Formula : C₁₇H₂₀N₄O₆
Molecular Weight : 376.36 g.mol⁻¹

3.2 Fluorescein

3.2.1 Introduction

In medicine, fluorescein is used in diagnostic techniques such as angiography, where it helps visualize blood flow in the retina.

In industry, it can be used to detect and locate leaks in products or packaging, detect leakage problems or trace underground watercourses.

3.2.2 Use

The microtube contains approx. 0.8 g Riboflavin, dose for one liter.

To observe fluorescence and create a measuring range, weigh and dilute the entire microtube in 1L of distilled water.


This solution will be used to study the compound's fluorescence response.

You can also carry out an assay session, using the stock solution to create a standard range.

Solution storage: Store for several days in a dark place.

Long-term storage will result in a drop in fluorescence intensity.

3.2.3 Results

Fluorescein Absorbance &

Fluorescence at 405nm
Fluorescence at 500 nm

Excitation

Absorbance maximum	450 nm – 510 nm
Spectrophotometer setting	405 nm – 500 nm

Emission

Fluorescence maximum	530 nm – 535 nm
----------------------	-----------------

3.2.4 Data

Name :	Fluorescein
CAS :	2321-07-5
Formula :	C ₂₀ H ₁₂ O ₅
Molecular Weight :	332.31 g.mol ⁻¹

3.3 Eosin

3.3.1 Introduction

In addition to its use in the medical field as a local antiseptic, eosin is a fluorescent dye often used in histology and cytology for tissue staining. It is particularly useful for highlighting cellular and tissue structures under the microscope. It is also used to diagnose leaks in networks and pipes, and to detect leaks in roofs and terraces..

3.3.2 Use

The microtube contains approx. 0.8 g Eosin, dose for one liter.

To observe fluorescence and set up a measuring range, accurately weigh and dilute the entire microtube in 1L of ethanol.



This solution will be used to study the fluorescence response of the compound.

An antiseptic solution dosing session (for example) can also be carried out, by creating a standard range from the stock solution.

Solution storage: Store for several days in a dark place.

Long-term storage will result in a drop in fluorescence intensity.

3.3.3 Results

Excitation

Absorbance maximum 525 nm
Spectrophotometer settings 500 nm

Emission

Fluorescence maximum 568 nm

3.3.4 Data

Name : Eosin Y (CI 45380)
CAS : 17372-87-1
Formula : $C_{20}H_6Br_4Na_2O_5$
Molecular Weight : 691.86 g.mol⁻¹

3.4 Rhodamine

3.4.1 Introduction

Rhodamine is a class of fluorescent dyes widely used in biology and chemistry to mark and visualize biological molecules and structures under a microscope. Also used to trace waterways, it is known for its lack of toxicity to flora and fauna..


3.4.2 Use

The bottle contains 0.2 g Rhodamine solution in 60 mL distilled water, dose for one liter. Fill the 60 mL vial to 1L with distilled water.

This solution will be used to study the fluorescence response of the compound. An antiseptic solution dosing session (for example) can also be carried out, by making a standard range from the stock solution.

Solution storage: store for several weeks/months in a dark place.

3.4.3 Results

Rhodamin Absorbance &

Fluorescence at 405nm

Excitation

Absorbance maximum 400-600 nm
Spectrophotometer settings 500 nm

Emission

Fluorescence maximum 612 nm

3.4.4 Data

Name : Rhodamin B
CAS : 81-88-9
Formula : $C_{28}H_{31}ClN_2O_3$
Molecular Weight : 479.01 g.mol⁻¹

3.5 Quinine

3.5.1 Introduction

In addition to its anti-malarial properties, quinine also has fluorescent properties. Extracted from the bark of the cinchona tree (a shrub native to South America), quinine is found in tonic beverages such as Schweppes® and Canada Dry®.


3.5.2 Application – Fluorescence

The microtube contains approx. 0.1 g Quinine, dose for one liter.
To observe fluorescence and set up a measuring range, accurately weigh and dilute the entire microtube in 1L of 1 M H₂SO₄.

This solution will be used to study the fluorescence response of the compound.
An energy drink dosing session can also be carried out, by making a standard range from the stock solution.

Solution storage: store for several weeks in a dark place.
Long-term storage will result in a drop in fluorescence intensity.

3.5.3 Results – Fluorescence

Excitation

Absorbance maximum 400 nm
Spectrophotometre settings 405 nm

Emission

Fluorescence maximum 667 nm

3.5.4 Use – Dosage


Preparation of a 100 mg.L⁻¹ solution (S₀)

- ⇒ Weigh the contents of the supplied microtube (approx. 0.1g)
- ⇒ Dilute in 1 L of 1 M H₂SO₄

Produce a range of Quinine standard solutions :

- ⇒ Concentrations to be achieved : from 10 to 100 mg.L⁻¹
- ⇒ Schweppes® contains a quinine concentration of around 70mg.L⁻¹, so there's no need to dilute it. However, be sure to degas the drink.

3.5.5 Results –Quinine determination in Schweppes®

The range enables the concentration of quinine in a spectrophotometer cell to be determined directly and graphically.

Here, we obtain a theoretical concentration of 71.6 mg.L⁻¹, in line with theory.

3.5.6 Data

Name :	Quinine hemisulfate monohydrate 98%
CAS :	61199-70-6
Formula :	C ₂₀ H ₂₄ N ₂ O ₂ .0,5H ₂ SO ₄ .H ₂ O
Molecular Weight :	391.47 g.mol ⁻¹

4. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

+33 (0)9 69 32 02 10