Observation

Stéréomicroscopes

Observation

Stereomicroscopes

**Réf :
571 265**

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0825 563 563 *

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
This high performance zoom stereomicroscopes is provided with a zoom system graduated from 1 to 4. The total magnification varies from x10 to x40.

These are meant to observe any biological or geological material, which cannot be observed with a microscope due to its thickness and/or its opacity.

2.Description

Caption :

- | | |
|---|----------------------------------|
| A Eyepiece tubes with eyeshades | E Illumination switches |
| B Dioptic adjustment | F Digital head |
| C Zoom adjustment | G Macrometric focusing |
| D LED incident and/or transmitted illumination | H Illumination adjustment |

3. Specifications

Head	Binocular
	45° inclined, 360° rotatable, with eyeshades
	1,3 Mp – USB – Windows 7 compatibility
Eyepieces	WF10x, Ø 30,5 mm
Interpupillar distance	55 - 75 mm
Dioptric adjustment	Both tubes
Objectives	Graduated zoom from 1 to 4 - Graduations : 1 - 1,5 - 2 - 3 - 4 with a separate knob
Magnification	X10 – x40
Field of view	23 mm
Working distance	91 mm
Stage	1 white/black and 1 transparent
Illumination	Both incident and transmitted LED with separate brightness adjustment
Power supply	4,5 V / 600 mA transformer
Dust cover	Yes

4. Unpacking and assembling

4.1. Unpacking the stereomicroscope

Carefully open the styrofoam packaging and check that all components have been taken out the packaging. Do not throw out the packaging before identifying all the components.

4.2. Assembling the microscope

When manipulating the components, and particularly optical parts, take care to not touch the lenses with your bare fingers or hands. Any fingerprints or grease spots will reduce the image quality.

After unpacking the support, install it on a stable working surface.

Make sure that the locking screws at the back are well fixed onto the metallic stand of the microscope. This operation is important because it prevents the microscope from accidentally sliding down along the stand.

Turn and fix the binocular body with the locking screw located on the left. Place the eyeshades onto the eyepieces.

Switch the instrument on by connecting the transformer to the electrical power supply. Before connecting the instrument and switching it on, make sure that the operating voltage of the stereomicroscope is the same as the power supply voltage.

4.3. Alignment and operation of the microscope

4.3.1. Interpupillary distance

Move the two eyepiece tubes until only one circular field is visible through the two eyepieces. If two circles appear, the interpupillary distance is too large, and if 2 overlapping circles appear, the interpupillary distance is too small.

4.3.2. Focusing the stereomicroscope

Choose the lowest magnification and adjust the image with the macrometric knob until having the clearest image. Increase the magnification step by step and readjust each time the image until having the required image.

4.3.3. Stage removing

In regards of what you will observe, you will prefer a clear or dark stage. With this stereomicroscope is provided 2 stages : the one in transparent frosted glass, the other is a reversible black and white stage. In order to remove this stage, unscrew without removing it the screw in front of the stage. Then, lift the stage in the back.

The transmitted illumination is only used with the frosted glass stage.

5. Digital camera install

5.1.

6. Cleaning and maintenance

Protection against dust : cover the instrument with the protective cover provided if it remains unused for a long period. Never leave the eyepiece tube open and unprotected. Either leave the eyepiece in the tube (which is recommended if the microscope is used very frequently) or cover it with wrapping paper or a protective cap. It is recommended that eyepieces and other optical accessories should be kept in a dry box to protect them from dust and humidity when they are not in use.

Protection against humidity : the instrument must be kept at a distance from every water source, pipe and sink. Humidity in the room in which the instrument is to be used must be as low as possible (the relative humidity must remain less than 70%). It is recommended that all optical accessories should be kept in a dry box when they are not being used. It is strongly recommended that a dehumidifier and air conditioning should be used continuously if the environment is too humid.

Cleaning : If there is any dust present on the optical surface, try to remove it with a fan or compressed air.

There are two recommended possible methods for removing fingerprints, grease or dust stains that cannot be removed by blowing with air :

- Blow gently onto the glass surface and then wipe with a clean cloth, paper for lenses or a cotton swab. Note that small cotton fibres may remain on the surface of the lens if a cotton swab is used.

- Use a cotton swab to clean the surface of the lens thoroughly. Do not use aggressive solvents.

The shelf life of the LED is of 5,000 hours. In case of non-function of the lighting and after checking that the batteries are charged, send the stereomicroscope back to Jeulin.

7. Technical support

The guarantee is 2 years. The equipment must be returned to our workshops.

If you need any repairs, adjustments or spare parts, please contact:

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