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## ENGLISH

### **Información importante**

All information, data ,comment and prompt in this manual is contain in normal scientific knowledge and common craftwork technology .

Please ready this manual book carefully before operator the Microtome; especially pay more attention to chapter Three of Safety direction and Warnings in other chapters .

Only eligible operator can manipulate instrument according to the direction of this manual. This instrument is just for tissue section of histology and pathology.

Only professional and eligible maintenance staff can maintain the microtome.

Protect and Maintain will be supply by our After-sale service department.

Because of the up to date of technical continuous, we have the right to make improvement of technical specifications without notification

Nameplate on back of instrument indicates the information of model, identification number and manufacture date

### Defintion

“Warning” ---- if violate operation sequence, there will be an serious accident and it will hurt the operator or instrument severe

“Attention”--- if violate operation sequence, maybe cause a error of operation and not work continually , even damage the instrument possible

## Index

|                                                                             |    |
|-----------------------------------------------------------------------------|----|
| 1. General description .....                                                | 25 |
| 2. Main technical parameters .....                                          | 25 |
| 3. Instructions for safe operation .....                                    | 26 |
| 3.1 Safety device of handwheel .....                                        | 26 |
| 3.2 Safety device of blade holder.....                                      | 27 |
| 4. Illustration of microtome .....                                          | 29 |
| 5. Unpacking and installation .....                                         | 30 |
| 5.1 Installing the microtome .....                                          | 30 |
| 5.2 Requirement of placement .....                                          | 30 |
| 5.3 Installing of waste tray .....                                          | 31 |
| 6. Orientation system of specimen .....                                     | 31 |
| 6.1 Embedding cassette clamp .....                                          | 32 |
| 6.2 Exchange between embedding cassette clamp and standard block clamp..... | 33 |
| 6.3 Installation of tissue specimen .....                                   | 33 |
| 7. Fix system of blade holder .....                                         | 34 |
| 7.1 Installation of the blade holder.....                                   | 34 |
| 7.2 Installation of the disposable blade .....                              | 35 |
| 7.3 Adjustment of cutting angle.....                                        | 36 |
| 7.4 Sway blade holder .....                                                 | 36 |
| 7.5 Safety precautions of blade holder operation .....                      | 37 |
| 8. Operation of slicing .....                                               | 37 |
| 8.1 Selection of slice thickness.....                                       | 37 |
| 8.2 Fast marching handwheel .....                                           | 38 |
| 8.3 Trimming .....                                                          | 39 |
| 8.4 Slicing .....                                                           | 40 |
| 8.5 Precautions of paraffin section .....                                   | 40 |
| 8.6 Change of specimen .....                                                | 41 |
| 8.7 Finish work.....                                                        | 41 |
| 9 Problems and solutions .....                                              | 42 |
| Problems.....                                                               | 42 |
| 10. Cleaning and maintenance .....                                          | 43 |
| 10.1 Disassembly of blade holder .....                                      | 43 |
| 10.2 Cleaning.....                                                          | 44 |
| 11. Storage and transportation .....                                        | 45 |

|                                               |    |
|-----------------------------------------------|----|
| 12. Packing list .....                        | 45 |
| 13. Maintenance and after sales service ..... | 46 |

## 1. General description

Advanced manual rotary Microtome is used widely for routine and laboratory slice test in medical colleges, Animal and Plant Inspection and Quarantine organization, Scientific research organization. It is used for slicing different level of thickness tissue specimen section for pathological diagnosis.

Advanced manual rotary microtome is latest products which be designed to adapt to the traditional operating habit of broad range of consumers. Such microtome absorbs the distillate of microtome technology of our company with history more than 10 years, combining with Structure character of mainstream microtome in domestic market, adoptting the advanced import roller guide shoe from Switzerland(SCHNEEBERGER), nimble feeding system without gap, and exempt lubrication, guarantee the high precision, increase of service life.

Advanced manual Rotary Microtome with novel structure, advanced features, excellent performance.

Precise and reliable Microfeed Mechanism with alarm system for stroke limit

The handwheel in the right can be locked in any position, and locked in highest position, provide the maximal security and convenient.

The handwheel in the lift can lead the specimen advance and retreat fastly and conveniently

Preparation forceps can be adjusted up/down left/right, trimming all kind shape of tissue specimen easily.

New type of blade holder can be moveing around and location easily, improving the work efficiency and Economic effects obviously by utilizing sharply cutting edge of the disposable blade

The section waste can be collected by wide tray at bottom, it's more environment-friendly

Blend into new design of ergonomics concept, bring the indefatigability and confort effect to operator

## 2. Main technical parameters

Environmental temperature +10°C — +40°C

Environmental moisture 35% — 80%

Section thickness range 0.5 — 60 μm

Slice thickness adjustment 0 — 2 μm increments 0.5 μm

2 — 10 μm increments 1 μm


|                                |                       |                            |
|--------------------------------|-----------------------|----------------------------|
|                                | 10 — 20 $\mu\text{m}$ | increments 2 $\mu\text{m}$ |
|                                | 20 — 60 $\mu\text{m}$ | increments 5 $\mu\text{m}$ |
| Horizontal specimen stroke     | 25 mm                 |                            |
| Vertical specimen stroke       | 59 mm                 |                            |
| Specimen orientation           |                       |                            |
| Horizontal                     | 8°                    |                            |
| Vertical                       | 8°                    |                            |
| Rotation                       | $\pm 90^\circ$        |                            |
| Movement range of blade holder |                       |                            |
| front and back                 | $\pm 25\text{mm}$     |                            |
| left and right                 | $\pm 20\text{mm}$     |                            |
| Height of cutting edge         | 105 mm                |                            |
| Dimensions (LxWxH)             | 460x276x78 mm         |                            |
| Net weight                     | 29 kg                 |                            |
| Gross weight                   | 40 kg                 |                            |

### 3. Instructions for safe operation

The safety devices on the instrument and accessory must not be modified or removed

#### 3.1 Safety device of handwheel

The spanner controls the handwheel, lock the handwheel at any position of rotation.

2. safety handwheel    2. rocking handle of handwheel

Before adjusting the section and specimen must be lock the handwheel

The handwheel can be locked in any direction with the locking spanner. The locking spanner can be set on the locked or released position on the right side of the base. Also can push the rocking handle at the top position to open or release the locked


Lock of handwhell

1. Convenient locking during manipulating——push the pry bar to the top, the handwheel can be locked at any position
2. Safety locked be used during the break——face the direction of the handle to push the hand wheel, shaking to the acme, the lockpin slid into lock tube aumomatic, then the handwheel been locked


### 3.2 Safety device of blade holder

Finger protection system which can flip up and down

warning: when operating or changing the specimen, knife and disposal blade, even during breaks, the cutting edge must be always covered with the knife guard



Before operating or changing the specimen or knife, and during breaks, it must always cover the cutting edge with the knife guard

Every knife carrier has different shape knife guard (1) it can be used to cover the cutting edge completely. Turn over the knife guard handle to cover the knife edge.

#### 4. Illustration of microtome

**1.** Waste tray **2.** Lock spanner of handwheel **3.** Handwheel **4.** The lock spanner of blade holder **5.** Knife carrier **6.** Lock spanner of blade **7.** Specimen clamp **8.** Scale indication **9.** Adjust knob of scale **10.** Blade holder **11.** Speed handwheel

## 5. Unpacking and installation

The equipment is sent in a wooden box with porexpan to avoid problems during the transport.


Take out the accessory box and the instruction manual.

Remove the around foam packing.

**Do not carry the instrument with the handle of the handwheel, or Blade Carrier device.**

### 5.1 Installing the microtome

Place the microtome on a stable laboratory table, and back away from the wall about 150mm.


### 5.2 Requirement of placement

1. The position must meet the following requirements:

2. Stable, no vibration laboratory table
3. No vibration floor
4. Room temperature always between +5°C~+40°C


### 5.3 Installing of waste tray

Insert the waste tray as the illustrated indicates.

## 6. Orientation system of specimen

The specimen orientation system consists of specimen clamp and angle adjusting device. The specimen clamp can be adjusted right and left or up and down to obtain the ideal slicing angle.


1. specimen clamp
2. screw for left and right
3. collet exchange handle
4. locking spanner of

gear block **5.** screw for up and down

Attention: The specimen clamping system has been installed and adjusted appropriately before leaving the factory. Users don't need to dismount and adjust by themselves.


## 6.1 Embedding cassette clamp

The instrument is furnished with a disposable embedding cassette clamp which is Unique design, it can rapid clamp the tissue specimen, improve the efficiency of work.it can be exchanged between embedding cassette clamp and wax block clamp, facing the requirement of different kinds of tissue specimen.

**1.** Gear block **2.** Moving jaw **3.** Regular jaw **4.** Spanner control open and close

## 6.2 Exchange between embedding cassette clamp and standard block clamp


1. Plastic embedding cassette 2. Embedding cassette universal clamp 3. Handle of clamp exchange 4. Standard block clamp

When exchange the different clamps, locked the handwheel, it is better to shift out the part of blade holder, then counterclockwise rotation and release handle of collet exchange, when the lock is loosened then you can extract the embedding cassette clamp or standard block clamp in the direction of the left. It can be reinstall in the reverse order after exchange, pay attention to insert wedge into the dovetail groove, then turn clamp and exchange spanner, just tightening the movable lock.

## 6.3 Installation of tissue specimen

1. When using the ordinary paraffin tissue block (1), release the knob of specimen clamp (3) to open the mouth of tongs (2), put in paraffin specimen and tighten it.

2. When using disposable plastic embedding cassettes, release open and close spanner of embedding cassette clamp (5), after putting in embedding cassette (4), turn back the open and close spanner for restoration


When exchanging specimen clamp or loading tissue specimen, the handwheel must be always locked and the cutting edge covered with the knife guard!

## 7. Fix system of blade holder

### 7.1 Installation of the blade holder

1. Pull out waste tray, release adjustment spanner of blade holder (1), push the dovetail groove of holder forward along with Rail that in the base of instrument. The locking block of blade holder (2) placed in the microtome base enters the T shape groove (4) of the blade holder.


2. Lock adjustment spanner (1), tighten the blade holder on the base of instrument and insert waste tray.

## 7.2 Installation of the disposable blade


1. Turn down the knife guard (2), release lock spanner of blade (1); at this point sickle clip will bounce automatically

2. Tuck the disposable blade from the side of direction (3) and lock spanner (1), turn back knife guard to the safety position.


### 7.3 Adjustment of cutting angle

1. Release tightening handle (1)
2. Turn the blade holder to adjust the angle of blade holder according to different specimen hardness and section thickness
3. Lock the tightening handle, tighten the blade holder on the base of blade holder

### 7.4 Sway blade holder

Blade holder(2) can move from side to side, have three anchor point: left, middle, right which take full advantage of blade to rough slice and fine cut effectively, Meanwhile it can save the loss of blade and improve the efficiency of work

## 7.5 Safety precautions of blade holder operation

Warning: when operating or changing the specimen, knife and disposal blade, even during breaks, the cutting edge must be always covered with the knife guard

- Be careful when adjusting the disposable blade, because the edge of blade is extremely sharp, inattention may cause injury to your hand
- Don't leave disposable blade and blade holder around, keep always unused disposable blade in the knife box.
- Blade will never allow to place upward
- Do not pick up in any case the falling blade and the blader holder with knife by hand

## 8. Operation of slicing

### 8.1 Selection of slice thickness

1. Choose the thickness by adjustment knob (1) on the top right corner of microtome. It can be rotated effectively to select each level for slice thickness

- Section thickness range: 0.5 — 60  $\mu\text{m}$

- Slice thickness adjustment: 0 — 2  $\mu\text{m}$


increments of 0.5  $\mu\text{m}$

2 — 10  $\mu\text{m}$

increments of 1  $\mu\text{m}$


|                       |                               |
|-----------------------|-------------------------------|
| 10 — 20 $\mu\text{m}$ | increments of 2 $\mu\text{m}$ |
| 20 — 60 $\mu\text{m}$ | increments of 5 $\mu\text{m}$ |

2. Window of slice thickness (2) display the selected slice thickness


## 8.2 Fast marching handwheel

1. The specimen can be rapidly approached by fast marching handwheel on the left side of microtome
2. Turning the Fast marching handwheel in the direction of arrow, specimen moves forwards.
3. Turning the Fast marching handwheel reverse, specimen moves backwards.
4. The turning will be tough when specimen moves to the limit position of stroke, if turn continues, it will feel slippery to turn the Fast marching handwheel

### 8.3 Trimming

1. Lock the handwheel safely
2. Clamp the specimen
3. Clamp disposable blade
4. Turn Fast marching handwheel to move specimen in the direction of knife edge or release the spanner of blade holder to move entire subassembly of blade holder allowing knife edge to get closer to specimen
5. Check all of gripping holder of specimen clamping device and blade holder, ensure clamp reliably
6. Select a slice thickness (such as 30um) or move small Fast marching handwheel tiny allowing the specimen and knife edge get closer and touch gradually until trim out a smooth surface of specimen by turnning the big Fast marching handwheel.

## 8.4 Slicing

1. Select the thickness by adjustment knob on the top right corner of microtome
2. Knife angle should be chosen starting from the smaller angle and increasing gradually until picking out a favorable slice angle
3. Turn the handwheel uniformly and lower the speed to cut hard specimen when slicing
4. Reticle mark on the right of blade holder is used to select different kinds of slice angle
5. Slice specimen by turning the handwheel uniformly

## 8.5 Precautions of paraffin section

Warning: before installing or exchanging specimen clamp and working clearance, handwheel must be locked and the cutting edge covered with the knife guard.

1. Lock handwheel

Warning: Clamp specimen prior to install blade

2. Install cool paraffin specimen on the specimen holder

Warning: pay more attention in case the knife or disposable blade is already installed; knife edge is very sharp, incorrect manipulation will cause a severe injury.

3. Turn Fast marching handwheel to move specimen to back position.
4. Insert the blade into blade holder from the side and clamp enough
5. Upturn the knife guard for finger protection, cover knife edge
6. Adjust the angle of slice (try it, range 0 to 3)
7. Allow the blade holder get closer to specimen as far as possible
8. Adjust the position of specimen as far as possible parallel with the edge
9. Release the handwheel

Warning: turn the handwheel uniformly when slicing. When slicing hard tissue specimen, lower the speed to turning handwheel in order to attain ideal slice.

10. Turn handwheel, start trimming
11. When attain the surface which you want to slice, stop trimming.
12. Select the desired thickness, turn the handwheel clockwise with a uniform speed and slice.

## **8.6 Change of specimen**

Warning: before installing or exchanging specimen clamp and working clearance, handwheel must be locked and the cutting edge covered with the knife guard.

1. Lock handwheel
2. Upturn the knife guard for finger protection, cover knife edge
3. Release specimen clamp and carry down specimen, change a new one

## **8.7 Finish work**

1. Lock handwheel
2. Release sickle clip and take out the blade from blade holder, put blade into knife box.
3. Carry down specimen from specimen clamp
4. Clean out all of slice waste
5. Clean instrument

## 9 Problems and solutions

| Problems                                                                                                                                                       | Possible causes                                                                                                                                                                                                                | Corrective action                                                                                                                                                                                                                                                                                                                                                                   |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Non-even of sections: The section thickness varies from one section to another. In extreme cases, sections are skipped meaning that a section is not obtained. | <ul style="list-style-type: none"> <li>-Insufficient knife/blade inclination; consequently the Clearance angle is too small.</li> <li>-Insufficient clamping of specimen and/or knife.</li> <li>-Blunt knife/blade.</li> </ul> | <ul style="list-style-type: none"> <li>-Systematically try several Clearance angle adjustments until the optimum angle is found.</li> <li>-Check if all clamps are locked and screws are tightened on the specimen and knife carrier systems. Retighten the clamps and screws if necessary.</li> <li>-Use a different part of the cutting edge or use a new knife/blade.</li> </ul> |
| Compressed sections<br>The sections are extremely compressed, wrinkled or jammed together.                                                                     | <ul style="list-style-type: none"> <li>-Blunt knife/blade.</li> <li>-Specimen too warm.</li> <li>-Clearance angle too big.</li> <li>-Sectioning speed too high.</li> </ul>                                                     | <ul style="list-style-type: none"> <li>- Use a different part of the cutting edge or use a new knife/blade.</li> <li>- Cool the specimen on a cold plate.</li> <li>- Clearance angle adjustment; systematically decrease the Clearance angle until the optimum adjustment is obtained.</li> <li>- Rotate the handwheel at a lower speed.</li> </ul>                                 |
| The knife "rings" on the cutting specimen when slicing hard specimens. Sections exhibit scratches and chatter marks.                                           | <ul style="list-style-type: none"> <li>-Sectioning speed too high.</li> <li>-Clearance angle too big.</li> </ul>                                                                                                               | <ul style="list-style-type: none"> <li>-Rotate the handwheel at a lower speed.</li> <li>-Clearance angle adjustment; systematically decrease the clearance angle until the optimum adjustment is obtained.</li> </ul>                                                                                                                                                               |


|                                                                      |                                                                                                                                                   |                                                                                                                                                                                                  |
|----------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                      | -Insufficient clamping of paraffin forceps or knife carrier.                                                                                      | -Check if all levers are locked and screws are tightened on the specimen and knife carrier systems. Retighten the clamps and screws if necessary.                                                |
| No more specimen advance and consequently no section produced.       | <ul style="list-style-type: none"> <li>- Specimen reached the front feed limit.</li> <li>- The coarse feed wheel cannot rotate freely.</li> </ul> | <ul style="list-style-type: none"> <li>- Turn the coarse feed wheel in the appropriate direction to move the specimen towards the clearance limit.</li> <li>- Remove the obstruction.</li> </ul> |
| Not lighting-up and no display when turn-on the power supply switch. | Fuse is burned                                                                                                                                    | Replace the fuse                                                                                                                                                                                 |

## 10. Cleaning and maintenance

Warning: lock handwheel before cleaning and maintenance instrument and unload blade holder and blade.

### 10.1 Disassembly of blade holder

1. Release the set lever of blade holder on microtome base and carry down subassembly of blade holder
2. Release set lever of blade (1), carry down blade carefully
3. Release set lever of blade (2), pull out blade holder in lateral direction
4. Release set lever of slide carriage (3), unload slide carriage of blade holder

## 10.2 Cleaning

1. Lock handwheel before cleaning the microtome
2. Remove slice waste
3. Unload waste tray, dump slice waste
4. Unload subassembly of blade holder, remove dirty paraffin from clamping part between disposable blade and blade holder carefully.
5. Cleaning the instrument just like Family routine Health, remove dirty paraffin around blade holder carefully.

6. Only use mild commercial detergents or soap solution for cleaning. Do not use solvent that contains acetone and benzene!

7. Use dry fabric or dipped by little detergent to clean the surface of the instrument.

8. Make sure, no detergent liquid enter the instrument during the cleaning. After cleaning, dry the surface of the instrument with a soft cloth.

9. Keep clean all parts of the instrument after use. Use the blade holder to move on the sliding rail, infuse the detergent lubricant. Good maintenance to the instrument will prolong the life of the instrument.

Warning: it is allowed to use normal commercial detergents or soap solution for disinfection. Do not use solvent that contains gasoline, acetone or benzene! It is not allowed liquid to flow inside the instrument

10. Once cleaning is over, install entire subassembly of blade holder and waste tray according to opposite procedure

11. Microtome can work long time without daily care, however, in order to operate normally within a long period, it is advisable to clean the instrument after end of work. Lubricate every month the parts of blade holder by lubricating oil provided with the equipment. Lock Handle shaft, slipper block, pressed compact dump rail etc, in order to reach the best result of use.

## 11. Storage and transportation

1. Environmental temperature:  $-5^{\circ}\text{C} \sim +55^{\circ}\text{C}$ ;

2. Environmental moisture:  $\leq 85\%$

3. Atmospheric pressure: 860hpa  $\sim$  1060hpa

4. Before transportation of this machine, please check if it's locked and if marks are complete. The microtome must be erectly transported!

## 12. Packing list

- Microtome Advanced
- Hexagon spanner
- User's manual
- Brush
- Waste tray
- Cassettes

- Metal clamp
- Handle replacement

### **13. Maintenance and after sales service**

If any damage, please do not carry out any repairs on your own, as this will invalidate the warranty. Qualified service engineers authorized by our company may only carry out repairs.

If the user operate strictly in accord with this manual, and obey the instruction, we will provide two-year warranty from the date of delivery. Within one warranty period, all damages caused by quality defects, we will be responsible for repairs or replacement. Besides, we will provide long time maintenance service.

Because of the up to date of technical continuous, we have the right to make improvement of technical specifications without notification