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System for the identification of *Enterobacteriaceae* in 4 hours

SUMMARY AND EXPLANATION

RapiD 20 E is a standardized system for the identification of *Enterobacteriaceae* in 4 hours, which uses 20 miniaturized biochemical tests chosen for their highly discriminant value and adapted to rapid interpretation. The complete list of those organisms that it is possible to identify with this system is given in the Identification Table at the end of this package insert.

PRINCIPLE

The RapiD 20 E strip consists of 20 microtubes containing dehydrated substrates. These tests are inoculated with a bacterial suspension which reconstitutes the media. During incubation, metabolism produces color changes that are either spontaneous or revealed by the addition of reagents.

The reactions are read according to the Reading Table and identification is obtained by referring to the Analytical Profile Index or using the identification software.

CONTENT OF THE KIT (Kit for 25 tests)

- 25 RapiD 20 E strips
- 25 incubation boxes
- 25 result sheets
- 1 package insert

COMPOSITION OF THE STRIP

The composition of the RapiD 20 E strip is given in the Reading Table of this package insert.

REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED

Reagents / Instrumentation

- API® NaCl 0.85 % Medium, 2 ml (Ref. 20 070)
- Reagents : VP 1 + VP 2 (Ref. 70 422)
JAMES (Ref. 70 542)
- Oxidase (Ref. 55 635*)
* reference not sold in certain countries : use an equivalent reagent.
- Mineral oil (Ref. 70 100)
- McFarland Standard (Ref. 70 900)
- RapiD 20 E Analytical Profile Index (Ref. 20 790) or **apiweb**™ identification software (Ref. 40 011) (consult bioMérieux)

Material

- Pipettes or PSIpettes
- Ampule protector
- Ampule rack
- General microbiology laboratory equipment

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI M29-A, *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline* - Current revision". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories - CDC/NIH - Latest edition", or to the regulations currently in use in each country.
- Do not use reagents past the expiry date.
- Before use, check that the packaging of the various components is intact.
- Do not use strips which have been damaged : cupules deformed, etc.
- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient history, the source of the specimen, colonial and microscopic morphology of the strain and, if necessary, the results of any other tests performed, particularly the antimicrobial susceptibility patterns.

STORAGE CONDITIONS

The strips should be stored at 2-8°C until the expiry date indicated on the packaging.

SPECIMENS (COLLECTION AND PREPARATION)

RapiD 20 E is not for use directly with clinical or other specimens.

The microorganisms to be identified must first be isolated on a suitable culture medium, preferably containing lactose (BCP, MacConkey, etc.), according to standard microbiological techniques.

INSTRUCTIONS FOR USE

Oxidase test

The oxidase test must be performed according to the manufacturer's instructions for use. The result should be recorded on the result sheet as it is an integral part of the final profile (21st identification test).

Selection of colonies

RapiD 20 E is to be used exclusively for members of the *Enterobacteriaceae*. In general, only oxidase-negative and Gram-negative bacilli should be tested with this strip.

NOTE : Certain Gram-negative bacilli not belonging to the *Enterobacteriaceae* and oxidase-positive (*Aeromonas* and *Vibrio*) are perfectly identifiable by RapiD 20 E. Clinical or bacteriological features will dictate whether to use this strip.

Preparation of the strip

- Prepare an incubation box (tray and lid). It is not necessary to add water to the tray.
- Record the strain reference on the elongated flap of the tray. (Do not record the reference on the lid as it may be misplaced during the procedure.)
- Remove the strip from its individual packaging.
- Place the strip in the incubation box.

Preparation of the inoculum

- Open an ampule of API® NaCl 0.85 % Medium (2 ml) as indicated in the paragraph "Warnings and Precautions" of the package insert for this product, or use any tube containing 2 ml of 0.85 % physiological saline without additives.
- Using a pipette or PSIpette, pick 1-4 well-isolated colonies of identical morphology from the agar plate either by suction or by successive touches. **Only young cultures should be used (18-24 hours old).**
- Carefully emulsify to achieve a homogeneous bacterial suspension with a turbidity equivalent to 0.5 McFarland. This suspension must be used immediately after preparation.
- Streak a purity plate with this suspension.

NOTE : In order for the RapiD 20 E tests to function properly, it is imperative that the inoculum density be adjusted to 0.5 McFarland. In particular, a weaker inoculum may lead to false negative results.

Inoculation of the strip

- With the same pipette, distribute the suspension into the tubes of the strip. To avoid the formation of bubbles at the base of the tubes, tilt the strip slightly forwards and place the tip of the pipette or PSIpette against the side of the cupule.
- For the CIT test, add 2 drops of the suspension (approx. 50 µl) to fill the tube and lower portion of the cupule.
- For the other tests, only fill the tubes (approx. 50 µl per tube). **The accuracy of the filling is very important.** Tubes insufficiently filled or overfilled may be the source of false positive and false negative results.
- For the underlined tests (LDC, ODC and URE), completely fill the cupule with mineral oil.
- Close the incubation box and incubate at 36°C ± 2°C for 4 - 4 ½ hours.

READING AND INTERPRETATION

Reading the strip

- After incubation, read the strip by referring to the Reading Table.
- Record all spontaneous reactions on the result sheet. For the fermentation tests (sugar substrates), a **green** color indicates the onset of acidification and should be considered as a **positive** reaction.
- Reveal the VP and IND tests by adding the corresponding reagents :
 - VP Test : add 1 drop of each of VP 1 and VP 2 reagents. Wait 5-10 minutes. A **red** color indicates a **positive** reaction to be recorded on the result sheet.
 - IND Test : add 1 drop of JAMES reagent. The reaction takes place immediately. A **pink** color indicates a **positive** reaction to be recorded on the result sheet.

NOTE : Do not put the lids back on the strips during the reading of the results.

Interpretation

Identification is obtained with the **numerical profile**.

- Determination of the numerical profile :
On the result sheet, the tests are separated into groups of 3 and a value 1, 2 or 4 is indicated for each. By adding together the values corresponding to positive reactions within each group, a 7-digit profile number is obtained. The oxidase reaction constitutes the 21st test and has a value of 4 if it is positive.
- Identification :
This is performed using the database (V3.1)
 - * with the Analytical Profile Index :
 - Look up the numerical profile in the list of profiles.
 - * with the **apiweb**™ identification software :
 - Enter the 7-digit numerical profile manually via the keyboard.

+	+	-	-	-	-	+	+	-	+	-	+	+	+	+	+	-	-	-	-	-
1	2	4	1	2	4	1	2	4	1	2	4	1	2	4	1	2	4	1	2	4
ONPG	LDC	ODC	URE	CIT	PPA	MNT	ESC	ARA	XYL	ADO	RHA	CEL	MEL	SAC	TRE	RAF	GLU	IND	VP	OX
3			0			4			5			6			7			1		

3 045 671 Escherichia coli

Further tests may be proposed in case of low discrimination. Refer to the identification software or Analytical Profile Index.

QUALITY CONTROL

The media, strips, and reagents are systematically quality controlled at various stages of their manufacture. For those users who wish to perform their own quality control tests with the strip, it is preferable to use the strain **1. *Escherichia coli* ATCC® 11775** or else one of the following strains :

2. *Proteus hauseri**** ATCC 13315 3. *Klebsiella pneumoniae* ssp *pneumoniae* ATCC 35657

ATCC : American Type Culture Collection, 10801 University Boulevard, Manassas, VA 20110-2209, USA.

	ONPG	LDC	ODC	URE	CIT	PPA	MNT	ESC	ARA	XYL	ADO	RHA	CEL	MEL	SAC	TRE	RAF	GLU	IND	VP
1.	+	+	+	-	-	-	-	-	+	-	-	+	-	+	-	+	-	+	+	-
2.	-	-	-	+	-	+	-	-	-	-*	-	-	-	-	+	-	-	+	+	-
3.	+	+	-	V	+	-	+	+	+	+	+	+	+	+	+	+	+	+	-	***

* This result may vary depending on the culture medium used.

** Weak positive reactions may occur.

*** Strain identified as *Proteus vulgaris* group with the RapiD 20 E strip.

Profiles obtained after culture of the strains on BCP agar (Bromocresol purple lactose agar).

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

LIMITATIONS OF THE METHOD

- The RapiD 20 E strip is intended for the identification of *Enterobacteriaceae*. A few other oxidase-positive and Gram-negative bacilli displaying both fermentative and oxidative metabolism (*Aeromonas*, *Vibrio*) give interpretable results. In general, only oxidase-negative and Gram-negative bacilli should be tested with this strip.
- The biochemical identification of *Salmonella*, *Shigella* as well as enteropathogenic *Escherichia coli* should be considered as presumptive and be confirmed by serology.
- The RapiD 20 E strip is intended for the identification of *Enterobacteriaceae* in 4 - 4 ½ hours. Respect the incubation time.
- A change or modification in the procedure indicated (inoculum too weak, tubes insufficiently filled or overfilled) may lead to false negative results, which could generate identification errors.
- Only pure cultures of a single organism should be used.

RANGE OF EXPECTED RESULTS

Consult the Identification Table at the end of this package insert for the range of expected results for the various biochemical reactions.

PERFORMANCE

2365 collection strains and strains of various origins belonging to species included in the database were tested :

- 95.52 % of the strains were correctly identified (with or without supplementary tests).
- 2.75 % of the strains were not identified.
- 1.73 % of the strains were misidentified.

WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

READING TABLE

TESTS	ACTIVE INGREDIENTS	QTY (mg/cup.)	REACTIONS/ENZYMES	RESULTS	
				NEGATIVE	POSITIVE
ONPG	2-nitrophenyl-βD-galactopyranoside	0.076	βD-galactosidase (Ortho-NitroPhenyl-βD-Galactopyranosidase)	colorless	pale yellow to bright yellow
<u>LDC</u>	L-lysine	0.46	Lysine DeCarboxylase	yellow-green to bluish-grey	blue to blue-violet
<u>ODC</u>	L-ornithine	0.28	Ornithine DeCarboxylase	yellow-green to bluish-grey	blue to blue-violet
<u>URE</u>	urea	0.32	UREase	yellow	pink / pink-violet
CIT	trisodium citrate	0.132	CITrate utilization	yellow to yellow-green	green to blue
PPA	4-nitrophenylalanine	0.024	Para-PhenylAlanine deaminase	colorless	orange-brown
MNT	sodium malonate	0.132	MaloNaTe utilization	yellow	green to blue
ESC	esculin ferric citrate	0.08 0.0236	β-glucosidase (ESCulin)	colorless	grey to black
ARA	L-arabinose	0.64	acidification (ARAbinose)	blue	green to yellow
XYL	D-xylose	0.6	acidification (XYLose)	blue	green to yellow
ADO	D-adonitol	0.64	acidification (ADOnitol)	blue	green to yellow
RHA	L-rhamnose	0.64	acidification (RHAmnose)	blue	green to yellow
CEL	D-cellobiose	0.64	acidification (CELlobiose)	blue	green to yellow
MEL	D-melibiose	0.64	acidification (MELibiose)	blue	green to yellow
SAC	D-saccharose (sucrose)	0.64	acidification (SACcharose)	blue	green to yellow
TRE	D-trehalose	0.64	acidification (TREhalose)	blue	green to yellow
RAF	D-raffinose	0.64	acidification (RAFfinose)	blue	green to yellow
GLU	D-glucose	0.64	acidification (GLUcose)	blue	green to yellow
IND	L-tryptophane	0.12	INDole formation	JAMES / immediate pale green-yellow pink	
VP	sodium pyruvate	0.08	acetoïn production (Voges-Proskauer)	VP 1 + VP 2 / 5-10 min colorless red	
OX	(see oxidase test package insert)	-	cytochrome-OXidase	(see oxidase test package insert)	

- The quantities indicated may be adjusted depending on the titer of the raw materials used.
- Certain cupules contain products of animal origin, notably peptones.

PROCEDURE
IDENTIFICATION TABLE
LITERATURE REFERENCES
INDEX OF SYMBOLS

p. I
p. II
p. III
p. IV

ATCC is a used, pending and/or registered trademark belonging to American Type Culture Collection.

bioMérieux® SA
au capital de 12 029 370 €
673 620 399 RCS LYON
69280 Marcy-l'Etoile / France
Tel. 33 (0)4 78 87 20 00
Fax 33 (0)4 78 87 20 90
<http://www.biomerieux.com>



bioMérieux, Inc
Box 15969,
Durham, NC 27704-0969 / USA
Tel. (1) 919 620 20 00
Fax (1) 919 620 22 11

Printed in France


bioMérieux, the blue logo, API, **apiweb** and RapiD 20 E are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

METHODOLOGIE / PROCEDURE / METHODIK / TECNICA / PROCEDIMENTO /
ΔΙΑΔΙΚΑΣΙΑ / METOD / METODE / METODYKA

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