

Benchcap

Metallic fixed work bench. Equipped with two height adjusting jacks

Meuble de support fixe en métal. Equipé de deux vérins pour le réglage vertical

Tisch aus Metall, höhenverstellbar


Mueble de soporte fijo en metal. Equipado con dos elevadores para ajuste vertical

Mobile di supporto fisso in acciaio elettrozincato con piedini regolabili


Product datasheet Fiche produit Ficha del producto Produktdatenblatt Scheda prodotto




Benchcap 320/321/391/481




320/321




391



320/321/391/481




481




	320	321	391	481
Structure / Structure / Struktur / Estructura / Struttura	Corrosion resistant electro-galvanized steel coated with anti-acid polymer Acier électro-zingué anti-corrosion gainé par enrobage de polymère thermo-durcissable anti-acide Elektroverzinkter Stahl, korrosionsfrei, durch thermisch härtbare Polymere umhüllt, säurefest Acero electro-galvanizado resistente a la corrosión recubierto de polímero termo-endurecedor anti-ácido Acciaio elettrozincato anticorrosione ricoperto con polimero termindurante antiacido			
Horizontal levelling / Mise à niveau / Einstellbarkeit / Nivelación / Aggiustabile	4 adjustable jacks / 4 vérins ajustables / 4 einstellbare Schraubenwinden / 4 gatos ajustables / 4 piedini aggiustabili			
Shelf (optional) / Étagère (Option) / Innenregal (Zubehör) / Estantería (opcional) / Ripiano (opzionale)	Sliding / Semi-extractible / Halbausziehbar / semiesträibile / Semi-estraibile			
Front support for installation of power sockets and taps Supports en façade pour installation prises électriques, vannes fluides et gaz Vorderleiste für die Installation von elektrischen Steckdosen, Ventilen für Flüssigkeiten und Gasen Soporte en fachada para instalación de tomas eléctricas, válvulas de fluido y gas Predisposizioni in facciata per installazione di prese elettriche, rubinetti acqua e gas	5	7	9	


Benchcap 392/483/633/714




392/483/633/714




392




483



633



714



	392	483	633	714
Structure / Structure / Struktur / Estructura / Struttura	Corrosion resistant electro-galvanized steel coated with anti-acid polymer Acier électro-zingué anti-corrosion gainé par enrobage de polymère thermo-durcissable anti-acide Elektroverzinkter Stahl, korrosionsfrei, durch thermisch härtbare Polymere umhüllt, säurefest Acero electro-galvanizado resistente a la corrosión recubierto de polímero termo-endurecedor anti-ácido Acciaio elettrozincato anticorrosione ricoperto con polimero termoindurante antiacido			
Horizontal levelling / Mise à niveau / Einstellbarkeit / Nivelación / Aggiustabile	4 adjustable jacks / 4 vérins ajustables / 4 einstellbare Schraubenwinden / 4 gatos ajustables / 4 piedini aggiustabili			
Shelf (optionnal) / Étagère (Option) / Innenregal (Zubehör) / Estanteria (opcional) / Ripiano (opzionale)	Sliding / Semi-extractible / Halbausziehbar / semiestraibile / Semi-estraiibile		Fixed	
Front support for installation of power sockets and taps Supports en façade pour installation prises électriques, vannes fluides et gaz Vorderleiste für die Installation von elektrischen Steckdosen, Ventilen für Flüssigkeiten und Gasen Soporte en fachada para instalación de tomas electricas, válvulas de fluido y gas Predisposizioni in facciata per installazione di prese elettriche, rubinetti acqua e gas	7	9	11	13