

Milieu SP1 et SP2 en poudre pour 1000 mL Réf. 114004 et 114010

Composition du milieu pour le Saintpaulia

Constituants	SP1 (mg/l) réf. 114004	SP2 (mg/l) réf. 114010
KNO₃	1900	950
NH₄NO₃	400	200
MgSO₄·7H₂O	370	185
CaCl₂·2H₂O	440	220
KH₂PO₄	170	85
MnSO₄·4H₂O	22,3	11,15
ZnSO₄·7H₂O	8,6	4,3
H₃BO₃	6,2	3,1
KI	0,83	0,415
Na₂MoO₄·2H₂O	0,25	0,125
CuSO₄·5H₂O	0,025	0,0125
CoCl₂·6H₂O	0,03	0,015
Fe, Na EDTA,	36,7	18,35
Inositol	100	100
thiamine HCl	0,5	0,5
ac. nicotinique	0,5	0,5
Pyridoxine, HCl	0,5	0,5
Ca panthoténate	0,5	0,5
biotine	0,01	0,01
A.I.A. (auxine)	0,2	
B.A.P. (cytokinine)	0,2	
SUCRE	20000	20000
GÉLOSE	7500	7150