

Mécanique

Ultrasons

Mechanics

Ultrasounds

Réf :
222 096

Français – p 1

English – p

Version : 2206

Récepteur Moduson®3

Moduson®3 Receiver

1. Description



Ce récepteur permet de recevoir les ultrasons d'un émetteur télémètre Moduson® 3 réf. 222 094 (ou d'un émetteur de la gamme Moduson® 2, mais nécessite tout de même la présence de l'émetteur télémètre Moduson®3 pour la liaison électrique).
Les ultrasons reçus sont transformés en signal électrique mesurable.
Sa tête renforcée et étanche en fait un outil pédagogique robuste et polyvalent.

2. Caractéristiques techniques

Sensibilité : ≥ 200 mVcc à 1.0 m de l'émetteur
Fréquence de travail : 40 ± 0.5 kHz
Distance de détection : 0.3 à 3.5 m
Directivité : 60°
Sortie : connecteur Jack 3.5mm
Connexion sur le boîtier de l'émetteur télémètre Moduson®3 réf.222 094
Dimensions boîtier : 57 x 44 x 30 mm

3. Service après-vente

La garantie est de 2 ans.

Pour toutes questions, contacter le **Support Technique** au **09 69 32 02 10***.

Si le matériel doit être retourné dans nos ateliers, pour toutes réparations ou demande de pièces détachées, veuillez contacter :

JEULIN – S.A.V.
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX France

09 69 32 02 10*

** Appel non surtaxé, prix d'un appel local*

1. Description



This receiver can receive ultrasounds from ultrasound telemeter transmitter Moduson® 3 part.222 094 (or a rangefinder Moduson® 2 transmitter but have to be connected on the telemeter transmitter Moduson®).

The received ultrasound is transformed into measurable electrical signal.

The water proof reinforced head is a robust and multipurpose educational tool.

2. Specifications

Sensibility: ≥ 200 mVcc à 1.0 m from the transmitter

Nominal frequency: 40 ± 0.5 kHz

Detectable range: 0.3 à 3.5 m

Directivity: 60°

Input with Jack 3.5mm

Output on the telemeter transmitter Moduson®3 part. 222 094

Dimensions: 57 x 44 x 30 mm

3. After-sales service

The device is under a 2-year guarantee.

For any repairs, contact the **Technical Support** on **+33 (0)2 32 29 40 23**.

The equipment should be returned in our workshops and for any repairs or spares parts, please contact:

JEULIN – TECHNICAL SUPPORT

468 rue Jacques Monod

CS 21900

27019 EVREUX CEDEX FRANCE

+33 (0)2 32 29 40 23

